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Medical Oxygen Cylinder & Accessories sterilization and cleaning procedure

a. Oxygen Cylinder:

1. Use CGA approved oxygen equipment cleaner/equivalent solution on cylinder body.
2. 5.25% Chlorine bleach (Clorox regular liquid bleach/equivalent), mix 1 part with 4 parts of clean water. This solution can be recommended for cleaning oxygen concentrator before every use.

b. Oxygen Mask:

1. Apart the mask.
2. Soak all parts in warm (35 degree Celsius) soapy water for 20 minutes.
3. Rinse thoroughly in running water.
4. Soak in a 1-part white vinegar/equivalent solution to 3-parts water solution for 20 minutes.
5. Rinse thoroughly and let air dry.
6. Store in PVC plastic bag.
7. Replace oxygen mask in every month.
8. Clean/use new oxygen mask before providing new patient.

c. Oxygen Concentrator:

1. Wipe outside cabinet once per week with damp cloth using oxygen equipment cleaner or equivalent solution wipe outside cabinet before changing every patient. **N.B: Unplug the unit prior to wiping.**
2. 5.25% Chlorine bleach (Clorox regular liquid bleach/equivalent), mix 1 part with 4 parts of clean water. This solution can be recommended for cleaning oxygen concentrator before every use.
3. Filter: Replace the filter unit every 4 months.

d. Nasal Cannula:

1. Wipe the cannula tips every day with a soapy cloth and rinse.
2. Replace every month.
3. Clean/use new nasal cannula before providing new patient.

e. Oxygen Tube:

1. Wipe the outside of the tube with damp cloth every week.
2. Do not soak the long tube in water as it is very difficult to take all water out.
3. Replace oxygen tube in every 4 month.
4. Clean/use new oxygen tube before providing new patient.

Humidifier Bottle:

1. Use only when supply is greater than 4 liters/min.
2. Throw the water out when the water level reaches the refill line on the outside bottle.
3. Wash the bottle once every week/ use new bottle patient, following the cleaning steps in oxygen mask section.
4. Replace the bottle in every month.

**WHO Recommends:

1. 70% ethyl alcohol to disinfect small areas between uses, such as reusable dedicated equipment
2. Sodium hypochlorite at .05% (equivalent to 5000 ppm) for disinfecting areas.

*CGA (Compressed Gas association) Guidelines on Corona Virus contamination:

1. Medical Gas containers that are contaminated must be cleaned.
2. Personnel involved in handling must receive training and instruction of safe handling, use appropriate PPE.
3. Cleaning procedures must be performed in designated areas with appropriate equipment and facilities using safe and efficient waste disposal system.
4. Methods and materials used to clean medical gas containers must be approved by the supplier prior to use.

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