

MR Imaging System

ECHELON Smart

Sentinel Customer Support

TECHNICAL GUIDE

--- Installation ---



L1E-HM0106-05

CONFIDENTIALITY: This material including drawings contains proprietary information of Hitachi, Ltd. and it shall not be disclosed to the third party or reproduced without the prior written permission of Hitachi, Ltd.
© Hitachi, Ltd. 2017, 2019. All rights reserved.

Introduction

Trademarks and registered trademarks:

Company names and product names used in this manual are the trademarks and registered trademarks of each company.

Cautions on exportation:

When exporting this equipment, be sure to check the Foreign Exchange and Foreign Trade Control Law and the regulations related to export control in the United States of America, and perform the necessary procedures.

Contact Hitachi or an authorized representative if you have further inquiries.

Revision history:

First edition : January 2017
Fourth edition : February 2019
Fifth edition : March 2019

Contents

Introduction	iii
Contents	v
Chapter1 Summary	1
Chapter2 System structure	3
Chapter3 System requirements	5
Chapter4 Connection setting	7
4.1 Network structure	7
4.2 Broadband router setting.....	7
4.3 System setting.....	8
4.3.1 Preparing setting information	9
4.3.2 IP address registration in Network Connection	9
4.3.3 Proxy IP address setting.....	12
4.3.4 Internet Setting	17
4.4 Time setting of a system	25
4.5 Creation of user account.....	26
4.6 Client PC setting.....	27
4.6.1 Adobe® Flash® Player installation	27
Chapter5 Connection to system server	29
Chapter6 How to start Sentinel Customer Support.....	31
6.1 Language setting.....	31
6.2 Proxy server setting	35
6.3 How to start	38
Chapter7 Operation way	41
7.1 Way of basic operation.....	41
7.1.1 Navigation Pane.....	43
7.1.2 Content Pane.....	45
7.1.3 Action Pane.....	46
7.1.4 Refresh.....	48

7.1.5 Help	49
7.1.6 Exporting URL	50
7.2 Group view	54
7.2.1 Basic operation	54
7.2.2 Finding devices	56
7.2.3 Daily Alert Summary	62
7.2.4 Recent Activities	63
7.2.5 Alerts.....	64
7.2.6 Devices	72
7.2.7 Device Profile	75
7.2.8 Readings	78
7.2.9 Activities.....	82
7.2.10 ECHELON Smart monitoring items	83
7.2.11 Monitoring event log (system) errors.....	113
7.2.12 Monitor Service	116
7.2.13 Tool & Manual.....	117
7.3 Manage alerts	121
7.4 Run command	134
7.4.1 Basic operation	134
7.4.2 ECHELON Smart command.....	137
7.5 Add readings.....	189
7.6 File Management.....	192
7.6.1 Refresh device files.....	192
7.6.2 Transfer files	194
7.6.3 Manage system files.....	202
7.6.4 View transfer log.....	216
7.7 Connect to device	218
7.7.1 Before Start	219
7.7.2 How to start.....	221
7.7.3 VNC viewer	226
7.7.4 Personal information protection mode.....	229
7.7.5 Service Tool & Diagnostic.....	231
7.7.6 How to acquire and display image.....	232
7.7.7 How to get library	242
7.8 Update software	254
7.8.1 Displaying [Software update] Dialog	254
7.8.2 Service pack applicable schedule.....	257
7.8.3 Application result confirmation	261
7.9 Get actual configurations.....	262

7.10 View chart.....	266
7.11 Report.....	273
7.11.1 View report	273
7.11.2 Using content.....	279
7.11.3 Using search function	280
7.11.4 Report save	283
7.11.5 Report print.....	286
7.12 Personal Configuration.....	287
7.13 View event log (Only for administrator account)	299
7.14 View software version	301
7.15 View software license (Only for administrator account)	303
Chapter8 How to quit Sentinel Customer Support.....	305

Chapter1 Summary

Sentinel Customer Support makes system use smooth and assists for quick recovery in case of trouble. It connects customer's system and support center on the internet and provides various supports. It fulfills customer expectations with scrupulous service such as 24 hour monitoring, failure message, operation support and so on.

Chapter2 System structure

Structure of Sentinel Customer Support is as follows.

Customer system and our support center are connected via the internet, and the remote service is offered.




Figure 2-1 System structure

Service contents

- Sentinel Customer Support Server System always monitors the system via the internet.
- If Sentinel Customer Support Server System detects errors of System, person in charge gets mail notification.
- Accessing from Service Offices to Sentinel Customer Support Server System makes it possible to monitor System monitoring status.
- Accessing from Service Offices to Sentinel Customer Support Server System enables remote access to System.
- Web browsing becomes available from systems and customer convenience such as accessing to support information is improved.

Chapter3 System requirements

The following system requirements are needed for this system.

Table 3-1 System requirements

Customer side	System software
	Remote service option
	A contract with Internet service provider
	Internet access line : ADSL40M or higher
Support Center	PC (Windows 7)
	CPU: Pentium4 Processor (or compatible CPU) 2GHz or higher*
	RAM: 1GB or higher*
	Flash Player9 or higher
	Internet access line : ADSL40M or higher
	Cellular telephone (which PC can be as an alternative to) and which can receive Email.
	Microsoft .NET Framework Version 3.5 or higher

NOTE There is some possibility that a dialog like the following is displayed if PC, which does not fulfill the above system requirements, is used in Support Center. Use PC which fulfills the system requirements.



Figure 3-1 Error dialog

Chapter4 Connection setting

The following description explains each setting which is needed for Sentinel Customer Support use.

4.1 Network structure

The following is an example of network structure in customer site if an internet circuit is newly built. Refer to Chapter 4.2 for network setting. Refer to Chapter 4.3 for system setting.




Figure 4-1 Network structure

4.2 Broadband router setting

Set up DHCP server on the side of LAN.




Figure 4-2 DHCP server setting of Broadband Router

4.3 System setting

The following explains system setting. If connection way is changed after that, report to a system administrator. Set up system once and then if connection way is changed, report to a system administrator.

Set up software according to the following operation flow.




Figure4-3 Operation flow

4.3.1 Preparing setting information.....

Inquire of customer LAN administrator for the following information in advance and prepare.

- Site information(Address and longitude latitude)
- IP address of system. **NOTE1**
- Subnet mask **NOTE1**
- Default gateway (IP address) **NOTE1**
- DNS server (IP address) **NOTE1**
- IP address and port number of HTTP proxy server (if proxy server is used for HTTP transfer on customer LAN.)
- User name and password of HTTP proxy server (if proxy server is used for HTTP transfer on customer LAN and the user is authenticated in such a proxy server.)

NOTE1 It is unnecessary for DHCP.

4.3.2 IP address registration in Network Connection

Perform network settings in Network Connection to connect system to LAN. (The following settings are unnecessary for DHCP.)

- (1) Way of starting setting tool
Log in system with service user and start [Network Connection] from [System Settings].




Figure4-4 System Settings menu

(2) System network setting

Configure the settings required for internet connection according to the following procedure.

- a) After finishing the operation that is described in “(1) starting set up tool”, confirm that the screen of figure 4-5 is displayed, select “Internet Protocol(TCP/IP)” (on the bottom of the list) from “This Connection use the following items” of “General” tab and click [Properties] button.




Figure4-5 Network Connections screen

- b) After a) operation, confirm that the screen of figure 4-6 is displayed, input the value which is specified by customer LAN administrator in IP address, Subnet mask, Default gateway, Preferred DNS server, Alternate DNS server and click [OK].




Figure4-6 TCP/IP Properties

(3) Quitting setting tool

Restart system after closing screen according to the following order.

4.3.3 Proxy IP address setting.....

Set up as follows if customer LAN connects to internet via proxy server. It is not necessary to set up if there is no proxy server on customer LAN.

- (1) How to starting setting tool
Start [Remote Service Setting] from [System Settings].




Figure4-7 System Settings menu

The following screen is displayed after Remote Service Setting boots.

The case where Proxy server setting is not configured.

The case where Proxy Server setting is configured. (Example.)

Figure4-8 Remote Service Setting screen

NOTE If Proxy Server setting is configured according to the procedure described in “(2) Proxy Server setting” on the next page, subsequently setting value at that time is indicated initially when Service Setting boots.

(2) Proxy Server setting

The following description explains input item.

Table4-1 Input item

Item name	Explanation
Proxy Server Address	Input IP address and port number of HTTP proxy server to go through at the time of accessing to internet Web server. It is not necessary to input if proxy server is not used for HTTP transfer on customer LAN. Make sure to input in the following format if inputting. Proxy server address : Port number (Example:192.168.0.199:8080)
User Name	Input the user name of the proxy server above. It is not necessary to input if the user is not authenticated in proxy server. Besides, it becomes impossible to input when Proxy Server Address is not input.
Password,Password(Verify)	Input password of the proxy server above. It is not necessary to input if the user is not authenticated in proxy server. Besides, it becomes impossible to input when Proxy Server Address is not input.
WS Server Address	Specify the IP address of <ws.sentinel-service.com> on VPN. Since this value is assigned automatically, input nothing in the usual environment.
RD Server Address	Specify the IP address of <rd.sentinel-service.com> on VPN. Since this value is assigned automatically, input nothing in the usual environment.
Static IP Route For Remote Service	When routing is required for the network connection to Sentinel servers, specify an IP address for the gateway or router to use when forwarding. Since routing settings are unnecessary in the usual environment, input nothing. And, this value cannot be inputted when the proxy server is specified.

Besides, input value needs to fulfill the following condition.

Table4-2 Input value condition

Item name	Input value condition
Proxy Server Address	Only alphabet, number, hyphen, colon, period (the character generally used as IP address or host name) can be used.
User Name	Only alphabet and number can be used. It is necessary to make sure to input when password is being input.
Password,Password(Verify)	Only alphabet and number can be used. It is necessary for same value to be input in both fields.

Input example.

Example1:

If proxy server is not used on customer LAN,
Proxy Server Address: Not input(**NOTE1**)

NOTE1 It becomes impossible to input other item if proxy server address is not input. In this case, program considers that all items are not input and executes processing.

Example2:

If proxy server whose IP address is 192.168.0.100 and whose port number is 8080 is used on customer LAN and the user is not authenticated in this proxy server,

Proxy Server Address: 192.168.0.100:8080

Other input item: Not input

Example3:

If proxy server whose host name is proxy01.XXX-hospital.com and whose port number is 80 is used on customer LAN and the user is authenticated in this proxy server,

Proxy Server Address: proxy01.XXX-hospital.com:80

UserName: The user name specified by customer LAN administrator.

Password,Password(Verify): The password specified by customer LAN administrator.

After completing input, click [Save] and save setting. After completely saving setting, the message that shows normal end or whether there is an error or not is displayed. If error message is displayed, review input value according to the following table.

Table4-3 Message list

Message content	Occurrence condition • Corrective action
Data saving is completed.	Occurrence condition: The case where setting has been completely and normally saved.
Character that can't be used for server name or IP address is included in the proxy server name.	Occurrence condition: Characters other than alphabet, figure, colon, period, and hyphen are being input in Proxy Server Address. Corrective action: Confirm input content of Proxy Server Address and input again.
Only alphabets and figures can use Proxy Server User Name.	Occurrence condition: Characters other than alphabet and figure are being input in User Name. Corrective action: Confirm input content of User Name and input again.
Only alphabets and figures can use Proxy Server Password.	Occurrence condition: Characters other than alphabet and figure are being input in Password. Corrective action: Input Password and Password (Verify) again.
When the Proxy Server Password is input, User Name is also necessary.	Occurrence condition: Password is being input but User Name is not being input. Corrective action: Review input content. (Input User Name or delete input value of Password and Password(Verify))
Value of the password and the verify password is different.	Occurrence condition: Password is in disagreement with Password(Verify) input value. (The case where either Password or Password(Verify) input is being input is also included.) Corrective action: Input Password or Password(Verify) again.

(3) Quitting Setting tool

Click [Close] or select [Exit] from [File] menu to quit this tool. In doing so, the following confirmation dialog is displayed.

The case where setting tool was quit without changing input value or after having saved setting.



Figure4-9 Remote Service Setting quitting confirmation popup1

If quitting: Click [Yes]

If continuing input: Click [No]

The case where setting tool was quit without saving setting.




Figure4-10 Remote Service Setting quitting confirmation popup2

If saving setting and quitting: Click [Yes]

If quitting without saving setting: Click [No]

The following dialog is displayed if setting is saved normally.




Figure4-11 Setting normal termination.


4.3.4 Internet Setting

Internet connectivity enables web browsing from MRI system and improves customers' convenience such as access to support information. Follow the next procedure and operate.


- 1 Click launcher [System Settings] button and select [Internet Setting] from menu.



- 2 "Internet Setting" starts up.




- 3 Input URL of our homepages in [Home page] [Address:], and input IP address of customer site proxy server in [Proxy server] [Address:], and input port number in [Port:] likewise.



NOTE1 If no proxy server, leave [Proxy server] blank.

NOTE2 If no port number is specified, leave [Proxy server] blank.


- 4 Click [OK] button.




NOTE Close Internet Browser once and confirm setting if setting was performed while Internet Browser was starting up.

After setting, start up an Internet Browser and confirm. Follow the next procedure and operate.

- 1 Click launcher [Help] button and choose [Internet Browser] from menu. Browser starts up.



- 2 If set homepage is displayed, setting is completed.



- 3 Like general web browsers, click  to close.

- Reference

The latitude longitude of the major city in the world is shown referring the latitude longitude input below.

Country	Capital	Latitude	Longitude
Afghanistan	Kabul	34.28North	69.11East
Albania	Tirane	41.18North	19.49East
Algeria	Algiers	36.42North	03.08East
American Samoa	Pago Pago	14.16South	170.43West
Andorra	Andorra la Vella	42.31North	01.32East
Angola	Luanda	08.50South	13.15East
Antigua and Barbuda	W. Indies	17.20North	61.48West
Argentina	Buenos Aires	36.30South	60.00West
Armenia	Yerevan	40.10North	44.31East
Aruba	Oranjestad	12.32North	70.02West
Australia	Canberra	35.15South	149.08East
Austria	Vienna	48.12North	16.22East
Azerbaijan	Baku	40.29North	49.56East
Bahamas	Nassau	25.05North	77.20West
Bahrain	Manama	26.10North	50.30East
Bangladesh	Dhaka	23.43North	90.26East
Barbados	Bridgetown	13.05North	59.30West
Belarus	Minsk	53.52North	27.30East
Belgium	Brussels	50.51North	04.21East
Belize	Belmopan	17.18North	88.30West
Benin	Porto-Novo (constitutional cotonou (seat of gvnt)	06.23North	02.42East
Bhutan	Thimphu	27.31North	89.45East
Bolivia	La Paz (adm.)/sucre (legislative)	16.20South	68.10West
Bosnia and Herzegovina	Sarajevo	43.52North	18.26East
Botswana	Gaborone	24.45South	25.57East
Brazil	Brasilia	15.47South	47.55West
British Virgin Islands	Road Town	18.27North	64.37West
Brunei Darussalam	Bandar Seri Begawan	04.52North	115.00East
Bulgaria	Sofia	42.45North	23.20East
Burkina Faso	Ouagadougou	12.15North	01.30West
Burundi	Bujumbura	03.16South	29.18East
Cambodia	Phnom Penh	11.33North	104.55East
Cameroon	Yaounde	03.50North	11.35East
Canada	Ottawa	45.27North	75.42West
Cape Verde	Praia	15.02North	23.34West
Cayman Islands	George Town	19.20North	81.24West
Central African Republic	Bangui	04.23North	18.35East
Chad	N'Djamena	12.10North	14.59East
Chile	Santiago	33.24South	70.40West

Country	Capital	Latitude	Longitude
China	Beijing	39.55North	116.20East
Colombia	Bogota	04.34North	74.00West
Comoros	Moroni	11.40South	43.16East
Congo	Brazzaville	04.09South	15.12East
Costa Rica	San Jose	09.55North	84.02West
Cote d'Ivoire	Yamoussoukro	06.49North	05.17West
Croatia	Zagreb	45.50North	15.58East
Cuba	Havana	23.08North	82.22West
Cyprus	Nicosia	35.10North	33.25East
Czech Republic	Prague	50.05North	14.22East
Democratic People's Republic of	P'yongyang	39.09North	125.30East
Democratic Republic of the Congo	Kinshasa	04.20South	15.15East
Denmark	Copenhagen	55.41North	12.34East
Djibouti	Djibouti	11.08North	42.20East
Dominica	Roseau	15.20North	61.24West
Dominica Republic	Santo Domingo	18.30North	69.59West
East Timor	Dili	08.29South	125.34East
Ecuador	Quito	00.15South	78.35West
Egypt	Cairo	30.01North	31.14East
El Salvador	San Salvador	13.40North	89.10West
Equatorial Guinea	Malabo	03.45North	08.50East
Eritrea	Asmara	15.19North	38.55East
Estonia	Tallinn	59.22North	24.48East
Ethiopia	Addis Ababa	09.02North	38.42East
Falkland Islands (Malvinas)	Stanley	51.40South	59.51West
Faroe Islands	Torshavn	62.05North	06.56West
Fiji	Suva	18.06South	178.30East
Finland	Helsinki	60.15North	25.03East
France	Paris	48.50North	02.20East
French Guiana	Cayenne	05.05North	52.18West
French Polynesia	Papeete	17.32South	149.34West
Gabon	Libreville	00.25North	09.26East
Gambia	Banjul	13.28North	16.40West
Georgia	T'bilisi	41.43North	44.50East
Germany	Berlin	52.30North	13.25East
Ghana	Accra	05.35North	00.06West
Greece	Athens	37.58North	23.46East
Greenland	Nuuk	64.10North	51.35West
Guadeloupe	Basse-Terre	16.00North	61.44West
Guatemala	Guatemala	14.40North	90.22West
Guernsey	St. Peter Port	49.26North	02.33West
Guinea	Conakry	09.29North	13.49West
Guinea-Bissau	Bissau	11.45North	15.45West
Guyana	Georgetown	06.50North	58.12West
Haiti	Port-au-Prince	18.40North	72.20West

Country	Capital	Latitude	Longitude
Heard Island and McDonald Islands		53.00South	74.00East
Honduras	Tegucigalpa	14.05North	87.14West
Hungary	Budapest	47.29North	19.05East
Iceland	Reykjavik	64.10North	21.57West
India	New Delhi	28.37North	77.13East
Indonesia	Jakarta	06.09South	106.49East
Iran (Islamic Republic of)	Tehran	35.44North	51.30East
Iraq	Baghdad	33.20North	44.30East
Ireland	Dublin	53.21North	06.15West
Israel	Jerusalem	31.71North	35.10West
Italy	Rome	41.54North	12.29East
Jamaica	Kingston	18.00North	76.50West
Jordan	Amman	31.57North	35.52East
Kazakhstan	Astana	51.10North	71.30East
Kenya	Nairobi	01.17South	36.48East
Kiribati	Tarawa	01.30North	173.00East
Kuwait	Kuwait	29.30North	48.00East
Kyrgyzstan	Bishkek	42.54North	74.46East
Lao People's Democratic Republic	Vientiane	17.58North	102.36East
Latvia	Riga	56.53North	24.08East
Lebanon	Beirut	33.53North	35.31East
Lesotho	Maseru	29.18South	27.30East
Liberia	Monrovia	06.18North	10.47West
Libyan Arab Jamahiriya	Tripoli	32.49North	13.07East
Liechtenstein	Vaduz	47.08North	09.31East
Lithuania	Vilnius	54.38North	25.19East
Luxembourg	Luxembourg	49.37North	06.09East
Macao, China	Macau	22.12North	113.33East
Madagascar	Antananarivo	18.55South	47.31East
Malawi	Lilongwe	14.00South	33.48East
Malaysia	Kuala Lumpur	03.09North	101.41East
Maldives	Male	04.00North	73.28East
Mali	Bamako	12.34North	07.55West
Malta	Valletta	35.54North	14.31East
Martinique	Fort-de-France	14.36North	61.02West
Mauritania	Nouakchott	20.10South	57.30East
Mayotte	Mamoudzou	12.48South	45.14East
Mexico	Mexico	19.20North	99.10West
Micronesia (Federated States of)	Palikir	06.55North	158.09East
Moldova, Republic of	Chisinau	47.02North	28.50East
Mozambique	Maputo	25.58South	32.32East
Myanmar	Yangon	16.45North	96.20East

Country	Capital	Latitude	Longitude
Namibia	Windhoek	22.35South	17.04East
Nepal	Kathmandu	27.45North	85.20East
Netherlands	Amsterdam/The Hague (seat of Gvnt)	52.23North	04.54East
Netherlands Antilles	Willemstad	12.05North	69.00West
New Caledonia	Noumea	22.17South	166.30East
New Zealand	Wellington	41.19South	174.46East
Nicaragua	Managua	12.06North	86.20West
Niger	Niamey	13.27North	02.06East
Nigeria	Abuja	09.05North	07.32East
Norfolk Island	Kingston	45.20South	168.43East
Northern Mariana Islands	Saipan	15.12North	145.45East
Norway	Oslo	59.55North	10.45East
Oman	Masqat	23.37North	58.36East
Pakistan	Islamabad	33.40North	73.10East
Palau	Koror	07.20North	134.28East
Panama	Panama	09.00North	79.25West
Papua New Guinea	Port Moresby	09.24South	147.08East
Paraguay	Asuncion	25.10South	57.30West
Peru	Lima	12.00South	77.00West
Philippines	Manila	14.40North	121.03East
Poland	Warsaw	52.13North	21.00East
Portugal	Lisbon	38.42North	09.10West
Puerto Rico	San Juan	18.28North	66.07West
Qatar	Doha	25.15North	51.35East
Republic of Korea	Seoul	37.31North	126.58East
Romania	Bucuresti	44.27North	26.10East
Russian Federation	Moskva	55.45North	37.35East
Rawanda	Kigali	01.59South	30.04East
Saint Kitts and Nevis	Basseterre	17.17North	62.43West
Saint Lucia	Castries	14.02North	60.58West
Saint Pierre and Miquelon	Saint-Pierre	46.46North	56.12West
Saint vincent and the Greenadines	Kingstown	13.10North	61.10West
Samoa	Apia	13.50South	171.50West
San Marino	San Marino	43.55North	12.30East
Sao Tome and Principe	Sao Tome	00.10North	06.39East
Saudi Arabia	Riyadh	24.41North	46.42East
Senegal	Dakar	14.34North	17.29West
Sierra Leone	Freetown	08.30North	13.17West
Slovakia	Bratislava	48.10North	17.07East
Slovenia	Ljubljana	46.04North	14.33East
Solomon Islands	Honiara	09.27South	159.57East
Somalia	Mogadishu	02.02North	45.25East

Country	Capital	Latitude	Longitude
South Africa	Pretoria (adm.) / Cap Town (Legislative) / Bloemfontein (Judicial)	25.44South	28.12East
Spain	Madrid	40.25North	03.45West
Sudan	Khartoum	15.31North	32.35East
Suriname	Paramaribo	05.50North	55.10West
Swaziland	Mbabane (Adm.)	26.18South	31.06East
Sweden	Stockholm	59.20North	18.03East
Switzerland	Bern	46.57North	07.28East
Syrian Arab Republic	Damascus	33.30North	36.18East
Tajikistan	Dushanbe	38.33North	68.48East
Thailand	Bangkok	13.45North	100.35East
The Former Yugoslav Republic of Macedonia	Skopje	42.01North	21.26East
Togo	Lome	06.09North	01.20East
Tonga	Nuku'alofa	21.10South	174.00West
Tunisia	Tunis	36.50North	10.11East
Turkey	Ankara	39.57North	32.54East
Turkmenistan	Ashgabat	38.00North	57.50East
Tuvalu	Funafuti	08.31South	179.13East
Uganda	Kampala	00.20North	32.30East
Ukraine	Kiev (Rus)	50.30North	30.28East
United Arab Emirates	Abu Dhabi	24.28North	54.22East
United Kingdom of Great Britain and Northern Ireland	London	51.36North	00.05West
United Republic of Tanzania	Dodoma	06.08South	35.45East
United States of America	Washington	47.30North	120.30West
United States of Virgin Islands	Charlotte Amalie	18.21North	64.56West
Uruguay	Montevideo	34.50South	56.11West
Uzbekistan	Tashkent	41.20North	69.10East
Vanuatu	Port-Vila	17.45South	168.18East
Venezuela	Caracas	10.30North	66.55West
Viet Nam	Hanoi	21.05North	105.55East
Yugoslavia	Belgrade	44.50North	20.37East
Zambia	Lusaka	15.28South	28.16East
Zimbabwe	Harare	17.43South	31.02East

4.4 Time setting of a system

Set up time of a system. After that, it synchronizes with NTP server automatically.




Figure4-12 Time setting

NOTE Only if time zone is properly set up and a margin of error is 1 hour or less, it can be synchronized with NTP server.

4.5 Creation of user account

Inform system administrator of the following items and request creating user account.

- Login user name
- Login user password
- Name
- Company
- Position
- User authority (for service maintenance/for reference)
- Department
- Time zone
- Email address (Cellular telephone)
- Email address (Office)
- Address
- Telephone number (Cellular telephone)
- Telephone number (Office)
- FAX number

4.6 Client PC setting

The following software is required for the client PC that performs remote service. Install according to the following procedure.


4.6.1 Adobe® Flash® Player installation

Connect to “<http://www.adobe.com/downloads/>” from browser, observe site screen, download and install a minimum of Flash Player 9 or later.


Flash Player is automatically installed if it is not installed when Sentinel customer support (<https://www.sentinel-service.com/qss/>) is connected to.

Observe the following procedure to install.


- 1 Log in to PC as Administrator authority user.
- 2 In case of Windows 7, the following display is on the upper side of browser if <https://www.sentinel-service.com/qss/> is specified in address.



- 3 Select and click [Install ActiveX Control].




- 4 The following dialog is displayed, so click [Yes].



- 5 Click [Install] on the following dialog.



- 6 Installation is completed and it is possible to connect to Sentinel Customer Support.



Chapter5 Connection to system server

After setting up system according to section 4.3, when connecting to internet line, server is connected to automatically and setting information is registered in server. Let system administrator know the following site information of system that has been newly connected and ask system administrator to register site information and confirm whether it has been possible to connect to server or not.

Table 5-1 Site information

System	Model name
	Serial number NOTE
Customer site	Site name
	Address
	Telephone number
	FAX number
	Email address
	Web site
Person in charge	Name of person in charge
	Position
	Organization name
	Address
	Telephone number
	FAX number
	Email address

NOTE The serial number can be confirmed by "Version Information" of "Help".



Figure5-1 Version Information


Chapter6 How to start Sentinel Customer Support

6.1 Language setting


It is necessary “English [en]” is selected as Web browser language, before Sentinel customer support starts.

Observe the following procedure to set up.


- 1 Open Internet Explorer and click [Tools]-[Internet Options].




- 2 Internet Options is displayed. Click [Languages].




- 3 Make sure that “English [en]” is indicated on the top of “Language”.




- 4 Click “Add” if “English [en]” is not indicated.




- 5 Search “English [en]” from “Add Language”.



- 6 Select “English [en]” and click “OK”.



- 7 Make sure that “English [en]” is indicated in “Language”, click “OK” to complete.




6.2 Proxy server setting


The following setting is needed if proxy server is used in a network of systems monitoring terminal PC that connects with Sentinel customer support.

NOTE It is unnecessary for DHCP.


- 1 Open Internet Explorer and click [Tools]-[Internet Options].




- 2 Internet Options is displayed. Click [Connections].




3 Click "LAN Setting".



4 The following dialog is displayed, so click "Advanced".



- 5 The following dialog is displayed, so tick "Use the same proxy server for all protocols".




- 6 Click "OK" to end.

6.3 How to start

- 1 Boot Internet Explorer.
- 2 Specify <https://www.sentinel-service.com/qss/> as address.


NOTE Registering this address to “Adding to Favorites” makes access easier from next time.

- 3 The following log-in screen is displayed.
Input user name and password, and click [Login].




- (1) User Name
Input a user name.
- (2) Password
Input a password.
- (3) Login
Input User Name, Password and click [Login].


- 4 If user name and password coincide, the main screen of the following Remote Service is displayed.



- 5 If user name and password do not coincide, input user name and password again.



- 6 "The session has expired." is indicated after a certain period of time without operation. Input user name and password again.



Chapter7 Operation way

7.1 Way of basic operation

Basic operation of Sentinel Customer Support is performed on internet screen. Operation is possible in the same way as usual internet screen.

Screen structure

- Screen is composed of right, left and middle parts.
- Panels switches display.
- Left and right parts of the screen display the functions for each application.
- If you click the items on the left and right sides of the screen, the content screen is switched to information, setup and operation screen.




Figure 7-1 Main screen


- (1) Navigation Pane
By switching panels, display and selection of [GROUPS] [DEVICE TYPES] are possible.
- (2) Action Pane
By switching panels, such operation as [SERVICE] and [CONFIGURATION] can be selected.
- (3) Content Pane
Detail operated with [Navigation Pane] and [Action Pane] is reflected.
- (4) Log in information.
User log-in information is indicated. Refresh, Sign out and Help operation is performed here.

7.1.1 Navigation Pane

By clicking  and  panels in Navigation Pane, [GROUPS] and [DEVICE TYPES] can be displayed.

[GROUPS] displays folder hierarchy of regionally divided groups. They are managed with folder hierarchy.

(Example)



NOTE To change group setting, administrator user authority is needed.




Figure7-2 GROUPS

[DEVICE TYPES] indicates the device type registered in system.

(Example) ECHELON Smart




Figure7-3 DEVICE TYPES

7.1.2 Content Pane

[Content Pane] indicates device summary and description about groups selected in [Navigation Pane].

Daily Alert Summary		Recent Activities			
Status	Count	Name	Service	Action	Date
Closed	230		ACM	poll	2007/06/14 18:10:05
Duplicate	459		MonitorManager	pushProperties	2007/06/14 18:05:04
Suppressed	44		MonitorManager	pushProperties	2007/06/14 18:05:04
			RemoteAccessManager	pushProperties	2007/06/14 18:05:04
			MonitorManager	pushProperties	2007/06/14 18:05:04
			MonitorManager	pushProperties	2007/06/14 18:05:04

Alerts ▾ Active			
Name	Status	Device	Date

Devices 9 Devices

Figure7-4 Content Pane

7.1.3 Action Pane

By clicking **SERVICE** and **CONFIGURATION** panels, [SERVICE] and [CONFIGURATION] items can be displayed in [Action Pane].

[SERVICE] indicates device monitoring functions and so on. [CONFIGURATION] indicates personal managing function and system configuration function.

NOTE To change [SYSTEM CONFIGURATION], administrator user authority is required.




Figure7-5 SERVICE




Figure7-6 CONFIGURATION

Table 7-1 Action Pane operation

Action Pane	List	Task	Description
SERVICE	ACTIONS	Manage Alerts	See Chapter7.3
		Manage Schedules	Not used.
		Run command	See Chapter 7.4.
		Add readings	See Chapter 7.5.
	TOTAL ACCESS	Connect to device	See Chapter 7.7.
	FILE MANAGEMENT	Refresh device files	See Chapter 7.6.1.
		Transfer files	See Chapter 7.6.2.
		Manage system files	See Chapter 7.6.3.
		View transfer log	See Chapter 7.6.4.
	SENTINEL UPDATE	Update software	See Chapter 7.8.
		Get actual configurations	See Chapter 7.9.
	REPORTS	View chart	See Chapter 7.10.
View report		See Chapter 7.11.	
CONFIGURATION	PERSONAL CONFIGURATION	Edit your profile	See Chapter7.12.
	SYSTEM CONFIGURATION	Manage devices	Administrator user authority is required.
		Manage customers	
		Manage alerts and alarms	
		Manage software updates	
		Manage connection types	
		Launch CMDDB editor	
	ADMINISTRATION	View TotalAccess	
		View event log	See Chapter 7.13.
		View software version	See Chapter 7.14.
View software license		See Chapter 7.15.	

7.1.4 Refresh

Refresh is screen reloading function.

Time interval to reload a screen can be chosen if [Refresh] button in upper-right corner is clicked. Choose “Now” if refreshing a screen immediately. Besides, choose “Every*****” if refreshing a screen automatically and periodically.





Figure7-7 Refresh

7.1.5 Help.....

It is possible to look up items and so on of the page opened through on-line help. If you click  on the upper right, on line help is displayed on another window.





Figure7-8 Help

7.1.6 Exporting URL

URL in the state where a group and a device are chosen can be exported, and [Groups] screen and [Devices] screen can be chosen from its address.






Figure7-9 Exporting URL

The following is a procedure for exporting URL.

- 1 Display a screen, ([Groups] screen or [Devices] screen), into which URL is desirable to export.



- 2 Click a device name (or group name) at the top of the content frame.




- 3 URL is exported.



- 4 Choose URL and copy the address from right-clicking menu.



- 5 Paste it to a text editor such as Notepad and save it.




The following is a procedure for logging in from exported URL.

1 Open the URL saved before by a web browser.

2 Log-in screen is displayed.

If this screen is added to “Add to Favorites”, log-in becomes convenient later.



- 3 If inputting a user name and password and logging in, a screen at the time of having been exporting URL is displayed.

The screenshot shows the Sentinel Customer Support web application. The browser is Microsoft Internet Explorer. The address bar displays <https://www.sentinel-service.com/ps/>. The page title is "Sentinel Customer Support". The user is signed in as "Hitachi".

The main content area displays "MRI_ECHELON_Smart_40925" with a "Matched configuration" table and a "Readings" table.

Matched configuration

Parameter	Value
Name	Type 5
Type 5	Type 5
V6.0D	V6.0D
DWI Package	DWI Package
BASG Imaging	BASG Imaging
Water Excitation	Water Excitation
RAPID	RAPID
CE-MRA Package	CE-MRA Package

Readings

Property	Value	Timestamp
Site/Image Count	110,357	Nov 29, 2016 4:23:05 PM
Site/Print Count	0	Nov 29, 2016 4:23:05 PM
Site/Scan Count	3,862	Nov 29, 2016 4:23:05 PM
Site/Patient Count	118	Nov 29, 2016 4:23:05 PM
Site/SNR(HMC)	366	Nov 23, 2016 12:16:43 PM
Site/SNR(HMC) Meas. Date	11/18/2016 22:23:50	Nov 23, 2016 12:16:43 PM
System/Time Zone	Tokyo Standard Time	Nov 29, 2016 4:23:05 PM
System/FE PC	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/Spectroscopy PC (Option)	11,003 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/PSC	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/ExtIO	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM

The right sidebar contains sections for SERVICE, ACTIONS, TOTAL ACCESS, FILE MANAGEMENT, SENTINEL UPDATE, REPORTS, and CONFIGURATION.

7.2 Group view

7.2.1 Basic operation

If you select any group from [GROUPS], list of the devices registered in its group and lower-level group is displayed and the latest operating condition and alert presence can be confirmed.






Figure7-10 Group selection

- (1) GROUPS
Any group can be selected from [GROUPS] and device registered in its group and lower-level group can be selected.
- (2) DEVICE TYPES
Any device type can be selected from [DEVICE TYPES] and its device type device can be selected.
- (3) Daily Alert Summary
By selecting group, Daily Alert Summary is displayed.
- (4) Recent Activities
By selecting group, [Recent Activities] is displayed.

- (5) Alerts
By clicking this window, [Alerts] screen is displayed.
- (6) Devices
By clicking this window, [Devices] screen is displayed.
- (7) Finding devices
By clicking , finding devices becomes possible.

7.2.2 Finding devices.....

Device Search enables finding devices.

By clicking  icon, Device Search screen is displayed.




Figure7-11 Finding devices

- (1) Name
Device search is performed with name.
Example) MRI_ECHELON_Smart_40925
- (2) Serial Number
Device search is performed with serial number.
Example) 40925
- (3) Advanced Search
This is used to set up for more advanced search.
- (4) Search
By clicking [Search] button, search begins.
- (5) Reset
By clicking [Reset] button, entry is deleted.
- (6) OK
By clicking [OK] button, selected device is indicated.
- (7) Cancel
By clicking [Cancel] button, this screen closes.

The following is device search procedure.

1 Fill in device name or serial number.

Name	MRI_ECHELON_Smart_40917
Serial Number	40917

If searching by system name, it is possible to search in a flexible way by using wildcard “*” or “%” to input a system name.

Table 7-2 Wild card description


To find...	Enter in Name field...
all devices whose names begin with the string "abc"	abc* or abc%
all devices whose names end with the string "abc"	*abc or %abc
all devices whose names include the string "abc"	*abc* or %abc%
only the device whose name is "abc"	abc

2 Click [Search] button.

3 Search result is listed on screen.

DEVICE SEARCH RESULT				1 Devices
Friendly Name	Name	Description	Serial Number	
MRI_ECHELON_...	MRI_ECHELON_...	Hitachi Hospital	40917	

- Select displayed device name, click [OK] button and its description is indicated.



The following is more advanced search procedure.

1 Click [Advanced Search] to display next screen.

Figure7-12 Advanced Search

- (1) Device type filter
Device search is performed with device type name.
- (2) Friendly Name
Friendly Name is performed with friendly name.
- (3) Device ID
Device search is performed with device ID.
- (4) Manufacturer **NOTE1**
Device search is performed with maker name.
- (5) Product **NOTE1**
Device search is performed with product name.
- (6) Model **NOTE1**
Device search is performed with model name.
- (7) Version **NOTE1**
Device search is performed with version name.
- (8) Account name
Device search is performed with account name.

- (9) Status
Device search is performed with status.
- (10) Description
Device search is performed with description.

NOTE1 If Device type filter is selected, manufacturer, product, model and version are automatically filled in.

2 Select or enter items to filter.

SEARCH CRITERIA

Search for Device Group User Configuration
(* = any string or character)

Name

Serial Number

+Advanced Search

Device type filter

Friendly Name

Device ID

Manufacturer

Product

Model

Version


Account name

Status

Description

UnAssigned Devices

3 Click [Search] button.



4 Search result is listed on screen.

DEVICE SEARCH RESULT			1 Devices
Friendly Name	Name	Description	Serial Number
MRI_ECHELON_...	MRI_ECHELON_...	Hitachi Hospital	40917

- 5 Select the device name being indicated, click [OK] button and its description is displayed.

The screenshot displays the Sentinel Customer Support web application in a Microsoft Internet Explorer browser. The page title is "Sentinel Customer Support" and the user is signed in as "Hitachi". The main content area shows the configuration for the device "MRI_ECHELON_Smart_40925".

Actual configuration information (Nov 29, 2016 4:25:44 PM):

Name	Value
Type 5	Type 5
V6.0D	V6.0D
DWI Package	DWI Package
BASG Imaging	BASG Imaging
Water Excitation	Water Excitation
RAPID	RAPID
CE-MRA Package	CE-MRA Package

Readings (299 Properties):

Reading Name	Value	Last Update
Site/Image Count	110,357	Nov 29, 2016 4:23:05 PM
Site/Print Count	0	Nov 29, 2016 4:23:05 PM
Site/Scan Count	3,862	Nov 29, 2016 4:23:05 PM
Site/Patient Count	118	Nov 29, 2016 4:23:05 PM
Site/SNR(HMC)	366	Nov 23, 2016 12:16:43 PM
Site/SNR(HMC) Meas. Date	11/18/2016 22:23:50	Nov 23, 2016 12:16:43 PM
System/Time Zone	Tokyo Standard Time	Nov 29, 2016 4:23:05 PM
System/FE PC	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/Spectroscopy PC (Option)	11,003 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/PSC	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM
System/ExtIO	0 [0:OK / Other:NG]	Nov 23, 2016 2:20:25 PM

The interface also includes a "GROUPS" sidebar with a list of states, a "SERVICE" sidebar with various management actions, and a "Readings" table with checkboxes for each reading. The bottom status bar shows "Last activity: FileRepositoryManager uploadFileNotify (Nov 29, 2016...)" and "0 Alerts".

7.2.3 Daily Alert Summary

By selecting any group from [GROUPS], [Daily Alert Summary] indicates the number of device alert registered in its group and lower-level group.




Figure7-13 Daily Alert Summary

- (1) Status
Alert status is indicated.
- (2) Count
The number of alert is indicated.

Table7-3 Alert status

Definitions of alert states		
Active	Pending	This is the alert request's initial state when it is created. The alert request is "pending" processing.
	Submitted	The processing of the alert request has begun; it has been submitted to the alert processor.
	Sent	The processing of the alert request is completed; it has been successfully sent to the alert processor.
Completed	Acknowledged	A user has reviewed the alert request and acknowledged that it exists.
	Closed	The alarm condition has been cleared or otherwise handled, and no more user action is required. The alert request is closed.
	Duplicate	The alert request pertains to the same alarm and alert as another alert request for the same device that is currently being processed. This alert request, therefore, is a duplicate and can be ignored.
	Suppressed	The alert escalator has exhausted all escalation levels for the alert request, and processing has reached the end. Therefore, all automatic processing is now suppressed.

7.2.4 Recent Activities

By selecting any group from [GROUPS], [Recent Activities] indicates operation detail of the device registered in its group and lower-level group.

Name	Service	Action	Date
MRI_ECHELON_Smar...	FileRepositoryManager	uploadFileNotify	Nov 29, 2016 4:28:4...
MRI_ECHELON_Smar...	RemoteAccess	executeAsync	Nov 14, 2016 11:24...

Figure7-14 Recent Activities

- (1) Name
Device name is indicated.
- (2) Service
Service content is indicated.
- (3) Action
Action content is indicated.
- (4) Date
Action date is indicated.

7.2.5 Alerts

By selecting any group from [GROUPS], [Alerts] indicates alert detail of device registered in its group and lower-level group.




Figure7-15 Alerts

- (1) Name
Alert name is indicated.
An icon on the side of Alert name indicates Alert status.
- | # | Icon | Description |
|---|---------------|--|
| 1 | (Orange icon) | It is displayed if Alert is being delivered. |
| 2 | (Yellow icon) | It is displayed if Alert is already delivered. |
| 3 | No icon | Alert is closed. Refer to Chapter 7.3. |
- (2) Status
Alert status name is indicated.
 - (3) Device
Device name is indicated.
 - (4) Date
Alert occurrence date is indicated.
 - (5) Display switching
[Table alerts] and [Map alerts] display are switched.
 - (6) The number of alert
The number of alert is indicated.
 - (7) Narrowing display down
By clicking ▼, indicating which status alert can be selected
Active status is indicated by default.

The following is alert detail status narrowing down procedure.

- 1 By clicking ▼ at the time of Table alerts, Alert status items are indicated as below.










Figure7-16 Alert status items to narrow down.

- 2 Select status to narrow down so that its status alert detail is indicated.

Alerts ▾ Active 1 Alert  

Name	Status	Device	Date
 Warning_Alert	Sent	MRI_ECHELON_...	2006/11/14 6:36:17



Alerts ▾ Suppressed 2393 Alerts  















Name	Status	Device	Date
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 7:58:44
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 6:58:46
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 6:58:45
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 5:13:38
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 5:13:37
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 5:13:36
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 3:30:42
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 3:30:31
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 2:09:29
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 1:45:50
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 0:44:45
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 0:44:43
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 0:44:42
 Warning_Alert	Suppressed	MRI_ECHELON_Smart...	2007/01/17 0:44:41


Figure7-17 Alert indication narrowing down

The following is mapping switch procedure.


- 1 Click  and [Map alerts] is on screen.


NOTE Map display is possible only if group is country or state unit.


- 2  means alert occurrence place.



- 3 Click  and alert occurrence device name is indicated.



- 4 Select  of indicated device and alert description is displayed.




NOTE Let administrator user know if alert coordinate is not accurately indicated and request for correcting coordinate.


By selecting device, Alerts occurring in device and Alarm are indicated.

Alerts: This is such processing as server sending e-mail to customers at the time of Alarm occurrence.

Alarms: Alarms is sent to server if the value exceeding threshold exists in monitoring property.

The following is Alarm indication procedure.



- 1 Select device to open [Alerts] window.
- 2 Click  so that [Alarms] is displayed instead and Alarms description is indicated.



Abstract	Description	Detail	Error code	Date
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 6:52:02
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/11/14 6:33:34
GCPreAmplifierSupplyFlu...	The amount of Fluxing G...	SL031000_GC_Pre_A...	ERROR	2006/11/14 6:33:34
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 5:44:29
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 5:06:28
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 4:16:31
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 3:54:44
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 2:20:49
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/13 20:10:33
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 5:27:10
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 5:11:52
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:54:20
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:29:35
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:04:24
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/11 3:26:03
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/11 3:08:57

Figure7-18 Alarms

- (1) Abstract
Alarms summary is indicated.
- (2) Description
Alarms description is indicated.
- (3) Detail
Alarms detail is indicated.
- (4) Error code
Error code is indicated.
- (5) Date
Alarm occurrence date is indicated.

- (6) Display narrowing down.
By clicking ▼, error code Alert to indicate can be selected. All are indicated by default.
 - (7) The number of Alarms
The number of Alarms is indicated.
 - (8) Display switch
If clicking  , [Alarms] and [Alerts] display are switched.
- 3 During [Alarms] display, by clicking ▼, error code is indicated as follows.




Figure7-19 Narrowing down Alarm error code item

- 4 If you select error code to narrow down, its error code alarm description is indicated.

The figure illustrates the process of narrowing down alarms by error code. The top screenshot shows a list of 101 alarms with columns for Abstract, Description, Detail, Error code, and Date. The bottom screenshot shows the same window filtered to show only 6 alarms with the error code 'Warning_Alarm'.

Abstract	Description	Detail	Error code	Date
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 6:52:02
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/11/14 6:33:34
GCPreAmpliferSupplyFlu...	The amount of Fluxing G...	SL031000_GC_Pre_A...	ERROR	2006/11/14 6:33:34
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 5:44:29
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 5:06:28
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 4:16:31
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 3:54:44
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/14 2:20:49
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/13 20:10:33
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 5:27:10
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 5:11:52
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:54:20
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:29:35
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/12 4:04:24
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/11 3:26:03
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/11 3:08:57
MPXH Power Supply +5...	Power supply voltage o...	SL140500_MPXH_Po...	ERROR	2006/11/11 2:33:10

Abstract	Description	Detail	Error code	Date
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/11/14 6:33:34
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/11/10 0:27:12
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/10/19 22:58:42
ShieldRoomTemperatur...	Temperature of shield r...	SL100000_Shield_Room...	WARNING	2006/10/19 22:58:42
PrimaryReturnFlux War...	The amount of Fluxing fi...	SL030700_Primary_Retu...	WARNING	2006/10/19 22:51:07
Scan PC HDD Temperat...	HDD temperature of the ...	Temperature_Celsius_S...	WARNING	2006/10/12 23:45:40

Figure7-20 [Alarms] display narrowing down

7.2.6 Devices

If you select any group from [GROUPS], [Devices] indicates the information about device registered under its group.





Figure7-21 Devices

- (1) Name
Device name is indicated.
- (2) Description
Site address description is indicated.
- (3) Device Type
Device type is indicated.
- (4) Online status
Device online status is indicated like the following table. It takes 3 minutes to refresh online status.

Table7-5 Device status

#	Status	Description
1	● (Red dot)	The device is offline.
2	⚡ (Lightning bolt)	Device is online with high-speed connection.
3	● (Green dot)	Device is online, but connection speed is slow.
4	● (Yellow dot)	The Enterprise cannot determine if the device is online or offline.
5	● (Gray dot)	The device is contactable.

- (5) Upload of files
The status of file upload is displayed. By clicking , uploaded “File Listing” can be referred to.


- (6) The number of devices
The number of devices registered in group is indicated.


If selecting any device, registered device summary, [Device Profile], [Readings], [Controls], [Activities] and [Alerts] can be confirmed.

The figure consists of two parts. The top part is a table titled 'Devices' with a '1 Devices' indicator in the top right corner. The table has three columns: 'Name', 'Description', and 'Device Type'. A red box highlights the 'Name' cell containing 'MRI_ECHELON Smart...'. Below the table is a large black arrow pointing downwards.

The bottom part is a screenshot of the 'Sentinel Customer Support' web application in Microsoft Internet Explorer. The browser address bar shows 'https://www.sentinel-service.com/tps/'. The page title is 'Sentinel Customer Support' and it indicates the user is signed in as 'Hitachi'. The main content area shows the selected device 'MRI_ECHELON_Smart_40925'. A 'Matched configuration' table compares 'Actual configuration information' (dated Nov 29, 2016 4:25:44 PM) with 'Value'. The configuration details include Name, Type 5, V6.0D, DWI Package, BASG Imaging, Water Excitation, RAPID, and CE-MRA Package. Below this is a 'Readings' section with a 'Select All' checkbox and a list of 299 properties, including Site/Image Count (110,357), Site/Print Count (0), Site/Scan Count (3,862), Site/Patient Count (118), Site/SNR(HMC) (366), and various system parameters. The bottom of the page shows 'Controls', 'Activities', and 'Alerts' sections.

Figure7-22 Device selection

By clicking , “File Listing” can be referred to.



All Device Files	Uploaded On	Last Modified
DICOM20061223132543.zip	2006/12/23 13:30:36	2006/12/23 13:31:36
User300032006122403330...	2006/12/24 15:37:16	2006/12/24 15:37:19
ServiceBackup.zip	2007/01/05 19:18:26	2007/01/05 19:22:52
diag_dump_info_MRI_ECH...	2007/01/17 20:03:03	2007/01/17 20:03:25
diag_dump_info_MRI_ECH...	2007/01/19 18:41:23	2007/01/19 18:41:44
diag_dump_info_MRI_ECH...	2007/01/22 18:01:08	2007/01/22 18:01:28

Figure7-23 File Listing

NOTE No files can be downloaded, etc from this screen. Refer to Chapter 7.6 for “File Management”.

7.2.7 Device Profile

If you select device, its device summary is indicated.






Figure7-24 Device summary

- (1) Device name
Device name is indicated.
- (2) Device information
Device information is indicated. Table indicates device online status. It takes 3 minutes to refresh online status.
- (3) Actual configuration information
The following software installed in device is indicated.
 - System software
 - Option
 - Service pack
 - Connection Type (Broad Band)
- (4) Image
Device image is displayed.
- (5) Device Profile
By clicking , Device Profile is indicated. The status of file upload is displayed. By clicking , uploaded “File Listing” can be referred to.
- (6) Assigned to
The person in charge of the device is displayed.

By clicking , [Device Profile] screen is displayed.




Figure7-25 Device Profile

(1) Device Profile

Table7-6 Device Profile items


1.	Name	Device name is indicated.
2.	Unique Device Name	Unique device name is indicated.
3.	Device Type	Device type is indicated.
4.	Device Time Zone	The time zone where the device is located.
5.	Description	Device location building name and so on is indicated.
6.	Address1	Device location place Address 1 is indicated.
7.	Address2	Device location place Address 2 is indicated.
8.	City	Device location place city is indicated.
9.	State or Province	Device location place state or province is indicated.

10.	Postal Code	Device location place postal code is indicated.
11.	Country	Device location place country is indicated.
12.	Unique Device ID	Device ID is indicated.
13.	Latitude	Device location place latitude is indicated.
14.	Longitude	Device location place longitude is indicated.
15.	Network address	Network address is indicated. NOTE: It is indicated only at the time of direct connection to server.
16.	Network accessible	Network access status is indicated. NOTE: It is indicated only at the time of direct connection to server.
17.	Serial Number	Serial number is indicated.
18.	Retry count	The number of Retry is indicated.
19.	Retry Interval	Retry interval is indicated.
20.	Current Poll interval	Current polling interval is indicated.
21.	Pending poll interval	The latest polling interval is indicated if configuration changes.

(2) Close

By clicking [Close] button, this screen closes.

(3) Help

By clicking , on-line help is displayed.

By clicking **Assigned to** ▼, the person in charge of the device is displayed.




Figure7-26 Assigned to

NOTE Please select the group that belongs the device when you select the device.

7.2.8 Readings

Device monitoring target value can be confirmed in [Readings] screen. Select device and [Readings] window opens.





Figure7-27 Readings

- (1) Property name
Property name is indicated. Tick at the time of individual reading.
- (2) Property graph
Graph is colored in blue if property value is normal value and it is colored in red at the time of aberrant value. No graph is displayed for the property having no value.
- (3) Value
Property value is indicated.
- (4) Date
Property value acquisition date is indicated.
- (5) Update
If  is clicked, the latest value of property whose is ticked is read.

By pointing the cursor to property, its description is indicated.




Figure7-28 Property selection


(1) Property	BE PC/HDD Free Space C:
(2) Value	68413829120 Byte
(3) Data Type	Integer
(4) Min	1073741824
(5) Max	74000000000

Figure7-29 Property example


- (1) Property
Property name is indicated.
- (2) Value
Property value and unit are indicated.
NOTE “Not available” is indicated at the time of no property value.
- (3) Data Type
Property value data type is indicated.
- (4) Min
Property minimum value is indicated.
- (5) Max
Property maximal value is indicated.

A value of device monitoring target is periodically acquired, but a value of the latest property can be acquired from Reading screen. Follow the next procedure to operate.

1 Confirm a device to read is online.




2 Tick of property whose value is acquired from Readings screen.




NOTE The number of to tick at once is 30 or less.

3 Click refresh button .



4 Put a mouse on Refresh button to indicate time until the latest value is acquired.



5 The latest value is acquired.

	73402363904	Byte	2007/01/23 9:08:22
	68424429568	Byte	2007/01/23 11:06:08
	1912	Hour	2007/01/23 9:08:22
	29	degC	2007/01/23 10:08:24

7.2.9 Activities

[Activities] screen indicates selected device activities.

Choose device and open [Activities] window.




Figure7-30 Activities

- (1) Service
Operated service name is indicated.
- (2) Action
Service action detail is indicated.
- (3) Date
Action date is indicated.
- (4) Last activity
The latest activity date is indicated.

7.2.10 ECHELON Smart monitoring items

Since ECHELON Smart, the following items are monitored. Notification mail is sent to customers if monitoring item value corresponds to warnings or error conditions. Furthermore, notification mail is also sent if errors occur to systems.

Table7-7 ECHELON Smart monitoring items

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Common					
Site/Image Count	-	-	-	-	1 hour
Site/Print Count	-	-	-	-	1 hour
Site/Scan Count	-	-	-	-	1 hour
Site/Patient Count	-	-	-	-	1 hour
Site/SNR(HMC)	-	-	-	-	24 hours
Site/SNR(HMC) Meas. Date	-	-	-	-	24 hours
System/Time Zone	-	-	-	-	24 hours
System/FE PC	-	0	≠0	-	1 hour
System/Spectroscopy PC (Option)	-	-	-	-	1 hour
System/PSC	-	0	≠0	-	1 hour
System/ExtIO	-	0	≠0	-	1 hour
System/HTTPS Response Time	Sec.	-	-	-	24 hours
MR-RFIP-1N					
System/TIP Board-1	-	0	≠0	-	1 hour
System/TIP Board-2	-	0	≠0	-	1 hour
System/TIP Board-3	-	0	≠0	-	1 hour
System/TIP Board-4	-	0	≠0	-	1 hour
Magnet/Magnet System(Alarm)	-	Normal	-	Fault (Quench Alarm, Shield Temp. Alarm, Helium Level Alarm, Pressure High Alarm or Pressure Low Alarm are other than Normal.)	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Magnet/Magnet System(Warning)	-	Normal	Warning (Helium Level Warning, Compressor Pressure Alarm, Compressor Helium Discharge Alarm, Compressor Helium After cooler Alarm or Compressor water outlet Alarm are other than Normal.)	-	1 hour
Magnet/Magnet System (Quench Alarm)	-	Normal	-	(Press ERDU)	1 hour
Magnet/Magnet System (Shield Temp. Alarm)	-	Normal	-	- (More than 60K.)	1 hour
Magnet/Magnet System (Helium Level Alarm)	-	Normal	-	- (Less than 56.7%)	1 hour
Magnet/Magnet System (Pressure High Alarm)	-	Normal	-	- (More than 4500Pa)	1 hour
Magnet/Magnet System (Pressure Low Alarm)	-	Normal	-	- (Less than 500Pa)	1 hour
Magnet/Gas-Bypass	-	-	-	-	1 hour
Magnet/Compressor System(Alarm)	-	Normal	Warning (Compressor Pressure Alarm, Compressor Helium Discharge Alarm, Compressor Helium After cooler Alarm or Compressor water outlet Alarm are other than Normal.)	-	1 hour
Magnet/Cable Disconnection (MR system - SVU)	-	Normal	-	-	1 hour
Magnet/Magnet Critical Error	-	Normal	-	Fault (In the case of no communication between SVU and MR system)	1 hour
Magnet/Magnet System (Helium Level Warning)	-	Normal	- (Less than 66.7%)	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Magnet/Compressor Pressure Alarm	-	Normal	- (Pressure: More than 2.55MPa Compressor control terminal #1, #2 more than 10MΩ)	-	1 hour
Magnet/Compressor Helium Discharge Alarm	-	Normal	- (Gass temp.: More than 90°C Compressor control terminal #3, #4 more than 10MΩ)	-	1 hour
Magnet/Compressor Helium After cooler Alarm	-	Normal	- (Gass temp.: More than 60°C Control terminal #5, #9 more than 10MΩ)	-	1 hour
Magnet/Compressor Water Outlet Alarm	-	Normal	- (Cooling water temp.: More than 60°C Compressor control terminal #10, #11 more than 10MΩ)	-	1 hour
Magnet/SVU (+DC48V Error)	-	Normal	- (Supply voltage abnormality in SVU)	-	1 hour
Magnet/SVU (+DC24V Error)	-	Normal	- (Supply voltage abnormality in SVU)	-	1 hour
Magnet/SVU (+DC5V Error)	-	Normal	- (Supply voltage abnormality in SVU)	-	1 hour
Magnet/SVU (-DC5V Error)	-	Normal	- (Supply voltage abnormality in SVU)	-	1 hour
Magnet/Li-Bat Low Alarm	-	Normal	-	- (Less than 11.2VDC)	1 hour
Magnet/Li-Bat High Alarm	-	Normal	- (More than 16.8VDC)	-	1 hour
Magnet/Supply Power of Helium Level	-	Normal	- (Supply voltage abnormality of helium measurement)	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Magnet/Measure Voltage of Helium Level	-	Normal	- (Supply voltage abnormality of helium measurement)	-	1 hour
Magnet/Circuit of pressure control	-	Normal	- (Pressure Heater error)	-	1 hour
Magnet/EIS Control	-	Normal	- (EIS error)	-	1 hour
Magnet/Helium Level	%	-	-	-	1 hour
Magnet/Pressure	Pa	-	-	-	1 hour
Magnet/Pressure Heater Output	mW	-	-	-	1 hour
Magnet/Bore Temp.	K	-	-	-	1 hour
Magnet/Shield Temp.	K	-	-	-	1 hour
Magnet/Li-Bat Voltage	VDC	-	-	-	1 hour
Magnet/SVU Version	-	-	-	-	1 hour
Magnet/SVU Serial Number(Upper)	-	-	-	-	1 hour
Magnet/SVU Serial Number(Lower)	-	-	-	-	1 hour
Magnet/Helium Warning Notification	-	Normal	- (Less than 68%)	-	1 hour
Sense/SENSE (Warning)	-	Normal	Warning	-	1 hour
Sense/SENSE (Fault)	-	Normal	-	Fault	1 hour
Sense/GC Flow Alarm	-	Normal	-	Fault	1 hour
Sense/GPA Flow Alarm	-	Normal	-	Fault	1 hour
Sense/Valve (Standby)	-	Normal	-	Standby	1 hour
Sense/Valve Abnormal	-	Normal	-	Fault	1 hour
Sense/Communication	-	Normal	-	Fault	1 hour
Sense/Compressor Flow Alarm	-	Normal	-	Fault	1 hour
Sense/Primary Supply Temp Alarm	-	Normal	-	Fault	1 hour
Sense/Primary Supply Temperature	degC	-	-	-	1 hour
Sense/GC Return Flow	l/min	-	-	-	1 hour
Sense/GPA Return Flow	l/min	-	-	-	1 hour
Sense/Compressor Return Flow	l/min	-	-	-	1 hour
Sense/Version	-	-	-	-	1 hour
Sense/PCB Power Supply(+3.3V)	V	-	-	-	1 hour
Sense/PCB Power Supply(+24V)	V	-	-	-	1 hour
Sense/PCB Temperature	degC	-	-	-	1 hour
Sense/DIPSW1-1	-	-	-	-	1 hour
Sense/DIPSW1-2	-	-	-	-	1 hour
Sense/DIPSW1-3	-	-	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Sense/DIPSW1-4	-	-	-	-	1 hour
Sense/DIPSW1-5	-	-	-	-	1 hour
Sense/DIPSW1-6	-	-	-	-	1 hour
Sense/DIPSW1-7	-	-	-	-	1 hour
Sense/DIPSW1-8	-	-	-	-	1 hour
Sense/DIPSW2-1	-	-	-	-	1 hour
Sense/DIPSW2-2	-	-	-	-	1 hour
Sense/DIPSW2-3	-	-	-	-	1 hour
Sense/DIPSW2-4	-	-	-	-	1 hour
Sense/DIPSW2-5	-	-	-	-	1 hour
Sense/DIPSW2-6	-	-	-	-	1 hour
Sense/DIPSW2-7	-	-	-	-	1 hour
Sense/DIPSW2-8	-	-	-	-	1 hour
Sense/Primary Supply Temp Warning	-	Normal	Warning	-	1 hour
Sense/Compressor Flow Warning	-	Normal	Warning	-	1 hour
Sense/GC Flow Warning	-	Normal	Warning	-	1 hour
Sense/GPA Flow Warning	-	Normal	Warning	-	1 hour
Sense/DIPSW3-1	-	-	-	-	1 hour
Sense/DIPSW3-2	-	-	-	-	1 hour
Sense/DIPSW3-3	-	-	-	-	1 hour
Sense/DIPSW3-4	-	-	-	-	1 hour
Sense/DIPSW3-5	-	-	-	-	1 hour
Sense/DIPSW3-6	-	-	-	-	1 hour
Sense/DIPSW3-7	-	-	-	-	1 hour
Sense/DIPSW3-8	-	-	-	-	1 hour
Sense/DIPSW4-1	-	-	-	-	1 hour
Sense/DIPSW4-2	-	-	-	-	1 hour
Sense/DIPSW4-3	-	-	-	-	1 hour
Sense/DIPSW4-4	-	-	-	-	1 hour
Sense/DIPSW4-5	-	-	-	-	1 hour
Sense/DIPSW4-6	-	-	-	-	1 hour
Sense/DIPSW4-7	-	-	-	-	1 hour
Sense/DIPSW4-8	-	-	-	-	1 hour
GC/GC Temperature	-	Normal	-	Fault	1 hour
RFPA/Power Supply	-	Normal	-	Fault	1 hour
RFPA/Peak Heat	-	Normal	-	Fault	1 hour
RFPA/Drive Power	-	Normal	-	Fault	1 hour
RFPA/Average Power Heat	-	Normal	-	Fault	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
RFPA/Peak Power	-	Normal	-	Fault	1 hour
RFPA/Average Power	-	Normal	-	Fault	1 hour
RFPA/Voltage Standing Wave Ratio	-	Normal	-	Fault	1 hour
RFPA/Gate Pulse Width	-	Normal	-	Fault	1 hour
RFPA/Gate Duty	-	Normal	-	Fault	1 hour
RFPA/RF Power Amplifier (Wait)	-	Normal	-	Wait	1 hour
RFPA/RF Power Amplifier (Standby)	-	Normal	-	Standby	1 hour
RFPA/Cable	-	Normal	-	Fault	1 hour
RFPA/Communication	-	Normal	-	Fault	1 hour
RFPA/Status	-	Normal	-	Fatal	1 hour
RFPA/Forward Power	W	0~25000	-	-	1 hour
RFPA/Reflected Power	W	0~25000	-	-	1 hour
RFPA/MEMORY BROKEN	-	Normal	-	Fault	1 hour
RFPA/OVER HEAT	-	Normal	-	Fault	1 hour
RFPA/DEVICE	-	Normal	-	Fault	1 hour
RFPA/HV LOW	-	Normal	-	Fault	1 hour
GPA/GPA Error	-	Normal	-	Fault	1 hour
GPA/X Axis Error	-	Normal	-	Fault	1 hour
GPA/Y Axis Error	-	Normal	-	Fault	1 hour
GPA/Z Axis Error	-	Normal	-	Fault	1 hour
GPA/Primary Power Error	-	Normal	-	Fault	1 hour
GPA/Communication PCB Error	-	Normal	-	Fault	1 hour
GPA/External Status Error	-	Normal	-	Fault	1 hour
GPA/External Clock Error	-	Normal	-	Fault	1 hour
GPA/Status Error	-	Normal	-	Fatal	1 hour
GPA/Communication Error	-	Normal	-	Fault	1 hour
GPA/Error Code Overflow	-	Normal	-	Fault	1 hour
GPA/GPA TYPE1	-	-	-	-	1 hour
GPA/GPA TYPE2	-	-	-	-	1 hour
GPA/Error Initial 1	-	-	-	-	1 hour
GPA/Error No 1	-	-	-	-	1 hour
GPA/Error Initial 2	-	-	-	-	1 hour
GPA/Error No 2	-	-	-	-	1 hour
GPA/Error Initial 3	-	-	-	-	1 hour
GPA/Error No 3	-	-	-	-	1 hour
GPA/Error Initial 4	-	-	-	-	1 hour
GPA/Error No 4	-	-	-	-	1 hour
GPA/Error Initial 5	-	-	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
GPA/Error No 5	-	-	-	-	1 hour
GPA/Error Initial 6	-	-	-	-	1 hour
GPA/Error No 6	-	-	-	-	1 hour
GPA/Error Initial 7	-	-	-	-	1 hour
GPA/Error No 7	-	-	-	-	1 hour
GPA/Error Initial 8	-	-	-	-	1 hour
GPA/Error No 8	-	-	-	-	1 hour
GPA/Error Initial 9	-	-	-	-	1 hour
GPA/Error No 9	-	-	-	-	1 hour
GPA/Error Initial 10	-	-	-	-	1 hour
GPA/Error No 10	-	-	-	-	1 hour
GPA/Error Initial 11	-	-	-	-	1 hour
GPA/Error No 11	-	-	-	-	1 hour
GPA/Error Initial 12	-	-	-	-	1 hour
GPA/Error No 12	-	-	-	-	1 hour
GPA/Error Initial 13	-	-	-	-	1 hour
GPA/Error No 13	-	-	-	-	1 hour
GPA/Error Initial 14	-	-	-	-	1 hour
GPA/Error No 14	-	-	-	-	1 hour
GPA/Error Initial 15	-	-	-	-	1 hour
GPA/Error No 15	-	-	-	-	1 hour
GPA/Error Initial 16	-	-	-	-	1 hour
GPA/Error No 16	-	-	-	-	1 hour
GPA/X Tune Parameter 1	-	-	-	-	1 hour
GPA/X Tune Parameter 2	-	-	-	-	1 hour
GPA/X Tune Parameter 3	-	-	-	-	1 hour
GPA/X Tune Parameter 4	-	-	-	-	1 hour
GPA/Y Tune Parameter 1	-	-	-	-	1 hour
GPA/Y Tune Parameter 2	-	-	-	-	1 hour
GPA/Y Tune Parameter 3	-	-	-	-	1 hour
GPA/Y Tune Parameter 4	-	-	-	-	1 hour
GPA/Z Tune Parameter 1	-	-	-	-	1 hour
GPA/Z Tune Parameter 2	-	-	-	-	1 hour
GPA/Z Tune Parameter 3	-	-	-	-	1 hour
GPA/Z Tune Parameter 4	-	-	-	-	1 hour
Temperature & Humidity/Shield Room Temperature	degC	5~32	<5 or >32	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Temperature & Humidity/Shield Room Humidity	%	0~80	>80	-	1 hour
Temperature & Humidity/Machine Room Temperature	degC	5~32	<5 or >32	-	1 hour
Temperature & Humidity/Machine Room Humidity	%	0~80	>80	-	1 hour
Temperature & Humidity/Shield Room Communication	-	Normal	Warning	-	1 hour
Temperature & Humidity/Machine Room Communication	-	Normal	Warning	-	1 hour
Decoupler/Decoupler	-	Normal	-	Fault	1 hour
KP/1KP	-	Normal	-	Fault	1 hour
KP/2KP	-	Normal	-	Fault	1 hour
KP/3KP	-	Normal	-	Fault	1 hour
KP/4KP	-	Normal	-	Fault	1 hour
Door/DOOR	-	CLOSE	-	-	1 hour
MPX/MRI System	-	Enable	-	-	1 hour
MPX/MR POWER	-	ON	-	-	1 hour
MPX/ECO MODE(GPA)	-	OFF	-	-	1 hour
MPX/ECO MODE(COMP)	-	OFF	-	-	1 hour
MPX/MPXH2 Power Supply(+24V)	V	22.6~26.4	<22.6 or >26.4	-	1 hour
MPX/MPXH2 Power Supply(+5V)	V	4.5~5.5	<4.5 or >5.5	-	1 hour
MPX/MPXH2 Power Supply(+3.3V)	V	2.97~3.63	<2.97 or >3.63	-	1 hour
MPX/MPXH2 Version	-	-	-	-	1 hour
MPX/DIPSW1-1	-	OFF	-	-	1 hour
MPX/DIPSW1-2	-	OFF	-	-	1 hour
MPX/DIPSW1-3	-	OFF	-	-	1 hour
MPX/DIPSW1-4	-	OFF	-	-	1 hour
MPX/DIPSW1-5	-	OFF	-	-	1 hour
MPX/DIPSW1-6	-	OFF	-	-	1 hour
MPX/DIPSW1-7	-	OFF	-	-	1 hour
MPX/DIPSW1-8	-	OFF	-	-	1 hour
MPX/DIPSW2-1	-	OFF	-	-	1 hour
MPX/DIPSW2-2	-	OFF	-	-	1 hour
MPX/DIPSW2-3	-	OFF	-	-	1 hour
MPX/DIPSW2-4	-	OFF	-	-	1 hour
MPX/DIPSW2-5	-	OFF	-	-	1 hour
MPX/DIPSW2-6	-	OFF	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
MPX/DIPSW2-7	-	OFF	-	-	1 hour
MPX/DIPSW2-8	-	OFF	-	-	1 hour
MPX/DIPSW3-1	-	OFF	-	-	1 hour
MPX/DIPSW3-2	-	OFF	-	-	1 hour
MPX/DIPSW3-3	-	OFF	-	-	1 hour
MPX/DIPSW3-4	-	OFF	-	-	1 hour
MPX/DIPSW3-5	-	OFF	-	-	1 hour
MPX/DIPSW3-6	-	OFF	-	-	1 hour
MPX/DIPSW3-7	-	OFF	-	-	1 hour
MPX/DIPSW3-8	-	OFF	-	-	1 hour
MPX/DIPSW16-1	-	OFF	-	-	1 hour
MPX/DIPSW16-2	-	OFF	-	-	1 hour
MPX/DIPSW16-3	-	OFF	-	-	1 hour
MPX/DIPSW16-4	-	OFF	-	-	1 hour
MPX/DIPSW16-5	-	OFF	-	-	1 hour
MPX/DIPSW16-6	-	OFF	-	-	1 hour
MPX/DIPSW16-7	-	OFF	-	-	1 hour
MPX/DIPSW16-8	-	OFF	-	-	1 hour
MPX/DIPSW17-1	-	OFF	-	-	1 hour
MPX/DIPSW17-2	-	OFF	-	-	1 hour
MPX/DIPSW17-3	-	OFF	-	-	1 hour
MPX/DIPSW17-4	-	OFF	-	-	1 hour
MPX/DIPSW17-5	-	OFF	-	-	1 hour
MPX/DIPSW17-6	-	OFF	-	-	1 hour
MPX/DIPSW17-7	-	OFF	-	-	1 hour
MPX/DIPSW17-8	-	OFF	-	-	1 hour
MR-IRCP-M1					
Magnet/Magnet System (Alarm)	-	Normal	-	Fault	1 hour
Magnet/Magnet System (Warning)	-	Normal	Warning	-	1 hour
Magnet/Magnet System (Quench Alarm)	-	Normal	-	-	1 hour
Magnet/Magnet System (Shield Temp. Alarm)	-	Normal	-	-	1 hour
Magnet/Magnet System (Helium Level Alarm)	-	Normal	-	-	1 hour
Magnet/Magnet System (Pressure High Alarm)	-	Normal	-	-	1 hour
Magnet/Magnet System (Pressure Low Alarm)	-	Normal	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Magnet/Gas-Bypass	-	-	-	-	1 hour
Magnet/Compressor System (Alarm)	-	Normal	Warning	-	1 hour
Magnet/Cable Disconnection (MR System - SVU)	-	Normal	-	-	1 hour
Magnet/Magnet Critical Error	-	Normal	-	Fault	1 hour
Magnet/Magnet System (Helium Level Warning)	-	Normal	-	-	1 hour
Magnet/Compressor Pressure Alarm	-	Normal	-	-	1 hour
Magnet/Compressor Helium Discharge Alarm	-	Normal	-	-	1 hour
Magnet/Compressor Helium After Cooler Alarm	-	Normal	-	-	1 hour
Magnet/Compressor Water Outlet Alarm	-	Normal	-	-	1 hour
Magnet/SVU (+DC48V Error)	-	Normal	-	-	1 hour
Magnet/SVU (+DC24V Error)	-	Normal	-	-	1 hour
Magnet/SVU (+DC5V Error)	-	Normal	-	-	1 hour
Magnet/SVU (-DC5V Error)	-	Normal	-	-	1 hour
Magnet/Li-Bat Low Alarm	-	Normal	-	-	1 hour
Magnet/Li-Bat High Alarm	-	Normal	-	-	1 hour
Magnet/Supply Power of Helium Level	-	Normal	-	-	1 hour
Magnet/Measure Voltage of Helium Level	-	Normal	-	-	1 hour
Magnet/EIS Control	-	Normal	-	-	1 hour
Magnet/Helium Level	%	-	-	-	1 hour
Magnet/Pressure	Pa	-	-	-	1 hour
Magnet/Pressure Heater Output	mW	-	-	-	1 hour
Magnet/Bore Temp.	K	-	-	-	1 hour
Magnet/Shield Temp.	K	-	-	-	1 hour
Magnet/Li-Bat Voltage	VDC	-	-	-	1 hour
Magnet/SVU Version	-	-	-	-	1 hour
Magnet/SVU_Serial_Number	-	-	-	-	1 hour
Magnet/Helium Warning Notification	-	Normal	-	-	1 hour
Magnet/Cold Head Warning Notification	-	Normal	Warning	-	1 hour
Magnet/Cold Head Error Notification	-	Normal	-	Fault	1 hour
Sense/SENSE (Warning)	-	Normal	Warning	-	1 hour
Sense/SENSE (Fault)	-	Normal	-	Fault	1 hour
Sense/Compressor Flow Alarm	-	Normal	-	Fault	1 hour
Sense/GC Flow Alarm	-	Normal	-	Fault	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Sense/GPA Flow Alarm	-	Normal	-	Fault	1 hour
Sense/Valve (Standby)	-	Normal	-	-	1 hour
Sense/Valve Abnormal	-	Normal	-	Fault	1 hour
Sense/Communication	-	Normal	-	Fault	1 hour
Sense/Primary Supply Temp Alarm	-	Normal	-	Fault	1 hour
Sense/Water Leak Alarm	-	Normal	-	Fault	1 hour
Sense/Primary Supply Temp Warning	-	Normal	Warning	-	1 hour
Sense/Compressor Flow Warning	-	Normal	Warning	-	1 hour
Sense/GC Flow Warning	-	Normal	Warning	-	1 hour
Sense/GPA Flow Warning	-	Normal	Warning	-	1 hour
Sense/Compressor Flow Connection	-	Normal	Warning	-	1 hour
Sense/GC Flow Connection	-	Normal	Warning	-	1 hour
Sense/GPA Flow Connection	-	Normal	Warning	-	1 hour
Sense/SENSE_Board_ID	-	-	-	-	24 hours
Sense/SENSE_Board_Work	-	-	-	-	24 hours
Sense/Board Revision	-	-	-	-	24 hours
Sense/Firmware Version	-	-	-	-	24 hours
Sense/Primary Supply Temperature	degC	-	-	-	1 hour
Sense/GC Return Flow	l/min	-	-	-	1 hour
Sense/GPA Return Flow	l/min	-	-	-	1 hour
Sense/Compressor Return Flow	l/min	-	-	-	1 hour
Sense/PCB Power Supply (+24V)	V	21.6~26.4	<21.6 or >26.4	-	1 hour
Sense/PCB Power Supply (+3.3V)	V	2.97~3.63	<2.97 or >3.63	-	1 hour
Sense/PCB Temperature	degC	5~32	<5 or >32	-	1 hour
Sense/DIPSW1-1	-	-	-	-	24 hours
Sense/DIPSW1-2	-	-	-	-	24 hours
Sense/DIPSW1-3	-	-	-	-	24 hours
Sense/DIPSW1-4	-	-	-	-	24 hours
Sense/DIPSW1-5	-	-	-	-	24 hours
Sense/DIPSW1-6	-	-	-	-	24 hours
Sense/DIPSW1-7	-	-	-	-	24 hours
Sense/DIPSW1-8	-	-	-	-	24 hours
Sense/DIPSW2-1	-	-	-	-	24 hours
Sense/DIPSW2-2	-	-	-	-	24 hours
Sense/DIPSW2-3	-	-	-	-	24 hours
Sense/DIPSW2-4	-	-	-	-	24 hours
Sense/DIPSW2-5	-	-	-	-	24 hours
Sense/DIPSW2-6	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
Sense/DIPSW2-7	-	-	-	-	24 hours
Sense/DIPSW2-8	-	-	-	-	24 hours
Sense/MIZU FW WDT Reset	-	OFF	-	-	1 hour
RFPA/Power Supply	-	Normal	-	Fault	1 hour
RFPA/Peak Heat	-	Normal	-	Fault	1 hour
RFPA/Drive Power	-	Normal	-	Fault	1 hour
RFPA/Average Power Heat	-	Normal	-	Fault	1 hour
RFPA/Peak Power	-	Normal	-	Fault	1 hour
RFPA/Average Power	-	Normal	-	Fault	1 hour
RFPA/Voltage Standing Wave Ratio	-	Normal	-	Fault	1 hour
RFPA/Gate Pulse Width	-	Normal	-	Fault	1 hour
RFPA/Gate Duty	-	Normal	-	Fault	1 hour
RFPA/RF Power Amplifier (Wait)	-	Normal	-	-	1 hour
RFPA/RF Power Amplifier (Standby)	-	Normal	-	-	1 hour
RFPA/Cable	-	Normal	-	Fault	1 hour
RFPA/Communication	-	Normal	-	Fault	1 hour
RFPA/Status	-	Normal	-	Fatal	1 hour
RFPA/Forward Power	W	0~25000	-	-	1 hour
RFPA/Reflected Power	W	0~25000	-	-	1 hour
RFPA/MEMORY BROKEN	-	Normal	-	Fault	1 hour
RFPA/OVER HEAT	-	Normal	-	Fault	1 hour
RFPA/DEVICE	-	Normal	-	Fault	1 hour
RFPA/HV LOW	-	Normal	-	Fault	1 hour
GPA/GPA Error	-	Normal	-	Fault	1 hour
GPA/X Axis Error	-	Normal	-	Fault	1 hour
GPA/Y Axis Error	-	Normal	-	Fault	1 hour
GPA/Z Axis Error	-	Normal	-	Fault	1 hour
GPA/Primary Power Error	-	Normal	-	Fault	1 hour
GPA/Communication PCB Error	-	Normal	-	Fault	1 hour
GPA/External Status Error	-	Normal	-	Fault	1 hour
GPA/Status Error	-	Normal	-	Fatal	1 hour
GPA/Communication Error	-	Normal	-	Fault	1 hour
GPA/Error Code Overflow	-	Normal	-	Fault	1 hour
GPA/HV ON Error	-	Normal	-	Fault	1 hour
GPA/External Clock Warning	-	Normal	-	Warning	1 hour
GPA/GPA_TYPE	-	-	-	-	24 hours
GPA/Error_Code1	-	-	-	-	1 hour
GPA/Error_Code2	-	-	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
GPA/Error_Code3	-	-	-	-	1 hour
GPA/Error_Code4	-	-	-	-	1 hour
GPA/Error_Code5	-	-	-	-	1 hour
GPA/Error_Code6	-	-	-	-	1 hour
GPA/Error_Code7	-	-	-	-	1 hour
GPA/Error_Code8	-	-	-	-	1 hour
GPA/Error_Code9	-	-	-	-	1 hour
GPA/Error_Code10	-	-	-	-	1 hour
GPA/Error_Code11	-	-	-	-	1 hour
GPA/Error_Code12	-	-	-	-	1 hour
GPA/Error_Code13	-	-	-	-	1 hour
GPA/Error_Code14	-	-	-	-	1 hour
GPA/Error_Code15	-	-	-	-	1 hour
GPA/Error_Code16	-	-	-	-	1 hour
GPA/X Tune Parameter 1	-	-	-	-	24 hours
GPA/X Tune Parameter 2	-	-	-	-	24 hours
GPA/X Tune Parameter 3	-	-	-	-	24 hours
GPA/X Tune Parameter 4	-	-	-	-	24 hours
GPA/Y Tune Parameter 1	-	-	-	-	24 hours
GPA/Y Tune Parameter 2	-	-	-	-	24 hours
GPA/Y Tune Parameter 3	-	-	-	-	24 hours
GPA/Y Tune Parameter 4	-	-	-	-	24 hours
GPA/Z Tune Parameter 1	-	-	-	-	24 hours
GPA/Z Tune Parameter 2	-	-	-	-	24 hours
GPA/Z Tune Parameter 3	-	-	-	-	24 hours
GPA/Z Tune Parameter 4	-	-	-	-	24 hours
KP/1KP	-	Normal	-	Fault	1 hour
KP/2KP	-	Normal	-	Fault	1 hour
KP/3KP	-	Normal	-	Fault	1 hour
KP/4KP	-	Normal	-	Fault	1 hour
Door/DOOR	-	CLOSE	-	-	1 hour
MRM/MRI System	-	Enable	-	-	1 hour
MRM/MR POWER	-	ON	-	-	1 hour
MRM/ECO MODE (COMP)	-	OFF	-	-	1 hour
MRM/MRM_Rack_ID	-	Normal	-	-	1 hour
MRM/MRM_Board_ID	-	Normal	-	-	1 hour
MRM/MRM_Board_Work	-	Normal	-	-	1 hour
MRM/Board Revision	-	Normal	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
MRM/FPGA Revision	-	Normal	-	-	1 hour
MRM/MRM FW Version	-	Normal	-	-	1 hour
MRM/SRM Module Link Error	-	Normal	-	-	1 hour
MRM/RFPA Unit Error	-	Normal	-	-	1 hour
MRM/GPA Unit Error	-	Normal	-	-	1 hour
MRM/SENSE Unit Error	-	Normal	-	-	1 hour
MRM/GANTRY Sensor Error	-	Normal	-	-	1 hour
MRM/MR ON	-	OFF	-	-	1 hour
MRM/MRM Firmware Startup Completion	-	OFF	-	-	1 hour
MRM/DIPSW1-1	-	OFF	-	-	24 hours
MRM/DIPSW1-2	-	OFF	-	-	24 hours
MRM/DIPSW1-3	-	OFF	-	-	24 hours
MRM/DIPSW1-4	-	OFF	-	-	24 hours
MRM/DIPSW1-5	-	OFF	-	-	24 hours
MRM/DIPSW1-6	-	OFF	-	-	24 hours
MRM/DIPSW1-7	-	OFF	-	-	24 hours
MRM/DIPSW1-8	-	OFF	-	-	24 hours
MRM/DIPSW2-1	-	OFF	-	-	24 hours
MRM/DIPSW2-2	-	OFF	-	-	24 hours
MRM/DIPSW2-3	-	OFF	-	-	24 hours
MRM/DIPSW2-4	-	OFF	-	-	24 hours
MRM/DIPSW2-5	-	OFF	-	-	24 hours
MRM/DIPSW2-6	-	OFF	-	-	24 hours
MRM/DIPSW2-7	-	OFF	-	-	24 hours
MRM/DIPSW2-8	-	OFF	-	-	24 hours
MRM/DIPSW3-1	-	OFF	-	-	24 hours
MRM/DIPSW3-2	-	OFF	-	-	24 hours
MRM/DIPSW3-3	-	OFF	-	-	24 hours
MRM/DIPSW3-4	-	OFF	-	-	24 hours
MRM/DIPSW3-5	-	OFF	-	-	24 hours
MRM/DIPSW3-6	-	OFF	-	-	24 hours
MRM/DIPSW3-7	-	OFF	-	-	24 hours
MRM/DIPSW3-8	-	OFF	-	-	24 hours
MRM/DIPSW4-1	-	OFF	-	-	24 hours
MRM/DIPSW4-2	-	OFF	-	-	24 hours
MRM/DIPSW4-3	-	OFF	-	-	24 hours
MRM/DIPSW4-4	-	OFF	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
MRM/DIPSW4-5	-	OFF	-	-	24 hours
MRM/DIPSW4-6	-	OFF	-	-	24 hours
MRM/DIPSW4-7	-	OFF	-	-	24 hours
MRM/DIPSW4-8	-	OFF	-	-	24 hours
MRM/DIPSW5-1	-	OFF	-	-	24 hours
MRM/DIPSW5-2	-	OFF	-	-	24 hours
MRM/DIPSW5-3	-	OFF	-	-	24 hours
MRM/DIPSW5-4	-	OFF	-	-	24 hours
MRM/DIPSW5-5	-	OFF	-	-	24 hours
MRM/DIPSW5-6	-	OFF	-	-	24 hours
MRM/DIPSW5-7	-	OFF	-	-	24 hours
MRM/DIPSW5-8	-	OFF	-	-	24 hours
MRM/DIPSW6-1	-	OFF	-	-	24 hours
MRM/DIPSW6-2	-	OFF	-	-	24 hours
MRM/DIPSW6-3	-	OFF	-	-	24 hours
MRM/DIPSW6-4	-	OFF	-	-	24 hours
MRM/DIPSW6-5	-	OFF	-	-	24 hours
MRM/DIPSW6-6	-	OFF	-	-	24 hours
MRM/DIPSW6-7	-	OFF	-	-	24 hours
MRM/DIPSW6-8	-	OFF	-	-	24 hours
MRM/DIPSW7-1	-	OFF	-	-	24 hours
MRM/DIPSW7-2	-	OFF	-	-	24 hours
MRM/DIPSW7-3	-	OFF	-	-	24 hours
MRM/DIPSW7-4	-	OFF	-	-	24 hours
MRM/DIPSW7-5	-	OFF	-	-	24 hours
MRM/DIPSW7-6	-	OFF	-	-	24 hours
MRM/DIPSW7-7	-	OFF	-	-	24 hours
MRM/DIPSW7-8	-	OFF	-	-	24 hours
MRM/DIPSW8-1	-	OFF	-	-	24 hours
MRM/DIPSW8-2	-	OFF	-	-	24 hours
MRM/DIPSW8-3	-	OFF	-	-	24 hours
MRM/DIPSW8-4	-	OFF	-	-	24 hours
MRM/DIPSW8-5	-	OFF	-	-	24 hours
MRM/DIPSW8-6	-	OFF	-	-	24 hours
MRM/DIPSW8-7	-	OFF	-	-	24 hours
MRM/DIPSW8-8	-	OFF	-	-	24 hours
MRM/Machine Room Temperature	degC	5~32	<5 or >32	-	1 hour
MRM/Machine Room Humidity	%	0~80	>80	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
MRM/MRM Power Supply (+12V)	V	10.8 ~ 13.2	<10.8 or >13.2	-	1 hour
MRM/MRM Power Supply (+5V)	V	4.5~5.5	<4.5 or >5.5	-	1 hour
MRM/MRM Power Supply (+3.3V)	V	2.97~3.63	<2.97 or >3.63	-	1 hour
MRM/MRM FW WDT Reset	-	OFF	-	-	1 hour
CONSOLE/CONSOLE_Board_ID	-	-	-	-	24 hours
CONSOLE/CONSOLE_Board_Work	-	-	-	-	24 hours
CONSOLE/Board Revision	-	-	-	-	24 hours
CONSOLE/FPGA Revision	-	-	-	-	24 hours
CONSOLE/Firmware Version	-	-	-	-	24 hours
CONSOLE/OPECALL	-	OFF	-	-	1 hour
CONSOLE/BUZZER	-	OFF	-	-	1 hour
CONSOLE/ECG	-	OFF	-	-	1 hour
CONSOLE/I2S Link Error	-	-	-	-	1 hour
CONSOLE/I2S Link Error (OPIF to EXIF)	-	Normal	-	Fault	1 hour
CONSOLE/I2S Link Error (EXIF to MRM)	-	Normal	-	Fault	1 hour
CONSOLE/I2S Link Error (MRM to SRM)	-	Normal	-	Fault	1 hour
CONSOLE/I2S Link Error (SRM to MRM)	-	Normal	-	Fault	1 hour
CONSOLE/I2S Link Error (MRM to EXIF)	-	Normal	-	Fault	1 hour
CONSOLE/I2S Link Error (EXIF to OPIF)	-	Normal	-	Fault	1 hour
CONSOLE/DIPSW1-1	-	-	-	-	24 hours
CONSOLE/DIPSW1-2	-	-	-	-	24 hours
CONSOLE/DIPSW1-3	-	-	-	-	24 hours
CONSOLE/DIPSW1-4	-	-	-	-	24 hours
CONSOLE/DIPSW1-5	-	-	-	-	24 hours
CONSOLE/DIPSW1-6	-	-	-	-	24 hours
CONSOLE/DIPSW1-7	-	-	-	-	24 hours
CONSOLE/DIPSW1-8	-	-	-	-	24 hours
CONSOLE/DIPSW2-1	-	-	-	-	24 hours
CONSOLE/DIPSW2-2	-	-	-	-	24 hours
CONSOLE/DIPSW2-3	-	-	-	-	24 hours
CONSOLE/DIPSW2-4	-	-	-	-	24 hours
CONSOLE/DIPSW2-5	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
CONSOLE/DIPSW2-6	-	-	-	-	24 hours
CONSOLE/DIPSW2-7	-	-	-	-	24 hours
CONSOLE/DIPSW2-8	-	-	-	-	24 hours
CONSOLE/Board_ID_(I.Box)	-	-	-	-	24 hours
CONSOLE/Board_Work_(I.Box)	-	-	-	-	24 hours
CONSOLE/I.Box_Serial_Number	-	-	-	-	24 hours
CONSOLE/I.Box Revision	-	-	-	-	24 hours
CONSOLE/OPIF FW WDT Reset	-	OFF	-	-	1 hour
SRM/General Status	-	-	-	-	1 hour
SRM/SRM CPU Status	-	-	-	-	1 hour
SRM/SRM Board Status	-	-	-	-	1 hour
SRM/PNLIF Board Detect Status	-	-	-	-	1 hour
SRM/TDRV Board Detect Status	-	-	-	-	1 hour
SRM/Gantry Monitor Power Status	-	-	-	-	1 hour
SRM/Manual Evacuate Sw Status	-	-	-	-	1 hour
SRM/Origin Position Detect Status	-	-	-	-	1 hour
SRM/Patient Fan Status	-	-	-	-	1 hour
SRM/Laser Status	-	-	-	-	1 hour
SRM/Illumination Status	-	-	-	-	1 hour
SRM/Release SW LED Status	-	-	-	-	1 hour
SRM/Remote Mode Transition Status	-	-	-	-	1 hour
SRM/CPU Error	-	Normal	-	Fault	1 hour
SRM/Limit Switch Sensor Error	-	Normal	-	Fault	1 hour
SRM/Gantry Panel Error	-	Normal	-	Fault	1 hour
SRM/PNLIF Communication Error	-	Normal	-	Fault	1 hour
SRM/TDRV Communication Error	-	Normal	-	Fault	1 hour
SRM/Up Down Moving Error	-	Normal	-	Fault	1 hour
SRM/Longitudinal Moving Error	-	Normal	-	Fault	1 hour
SRM/ExTmpHum Sensor Error	-	Normal	-	Fault	1 hour
SRM/SRM_Board_ID	-	-	-	-	24 hours
SRM/SRM_Board_Work	-	-	-	-	24 hours
SRM/Board Revision	-	-	-	-	24 hours
SRM/FPGA1 Revision	-	-	-	-	24 hours
SRM/FPGA2 Revision	-	-	-	-	24 hours
SRM/FW Version	-	-	-	-	24 hours
SRM/CAN I/F Version	-	-	-	-	24 hours
SRM/DSW2-1	-	-	-	-	24 hours
SRM/DSW2-2	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
SRM/DSW2-3	-	-	-	-	24 hours
SRM/DSW2-4	-	-	-	-	24 hours
SRM/DSW2-5	-	-	-	-	24 hours
SRM/DSW2-6	-	-	-	-	24 hours
SRM/DSW2-7	-	-	-	-	24 hours
SRM/DSW2-8	-	-	-	-	24 hours
SRM/DSW3-1	-	-	-	-	24 hours
SRM/DSW3-2	-	-	-	-	24 hours
SRM/DSW3-3	-	-	-	-	24 hours
SRM/DSW3-4	-	-	-	-	24 hours
SRM/DSW3-5	-	-	-	-	24 hours
SRM/DSW3-6	-	-	-	-	24 hours
SRM/DSW3-7	-	-	-	-	24 hours
SRM/DSW3-8	-	-	-	-	24 hours
SRM/DSW4-1	-	-	-	-	24 hours
SRM/DSW4-2	-	-	-	-	24 hours
SRM/DSW4-3	-	-	-	-	24 hours
SRM/DSW4-4	-	-	-	-	24 hours
SRM/SRM Power Supply (+24V)	V	21.6 ~ 26.4	<21.6 or >26.4	-	1 hour
SRM/SRM Power Supply (+12V)	V	10.8 ~ 13.2	<10.8 or >13.2	-	1 hour
SRM/SRM Power Supply (+5V)	V	4.5 ~ 5.5	<4.5 or >5.5	-	1 hour
SRM/SRM Power Supply (+3.3V1)	V	2.97 ~ 3.63	<2.97 or >3.63	-	1 hour
SRM/SRM Power Supply (+3.3V2)	V	2.97 ~ 3.63	<2.97 or >3.63	-	1 hour
SRM/Longitudinal Full Stroke Value	mm	-	-	-	24 hours
SRM/Longitudinal Absolute Position Value	mm	-	-	-	1 hour
SRM/TABLE_General_Status	-	-	-	-	1 hour
SRM/TABLE_UP/DOWN_Position	-	-	-	-	1 hour
SRM/TABLE_IN/OUT_Position	-	-	-	-	1 hour
SRM/OPECALL DETECT LEVEL	-	-	-	-	1 hour
PNLIF/PNLIF Board Detect Status	-	-	-	-	1 hour
PNLIF/PNLIF Serial Communication Error	-	-	-	-	1 hour
PNLIF/PANEL-R Detect Status	-	-	-	-	1 hour
PNLIF/PANEL-L Detect Status	-	-	-	-	1 hour
PNLIF/Check Sum Error (SRM to PNLIF)	-	-	-	-	1 hour
PNLIF/Com-ID Error (SRM to PNLIF)	-	-	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
PNLIF/Link Up Status (SRM to PNLIF)	-	-	-	-	1 hour
PNLIF/Check Sum Error (PNLIF to SRM)	-	-	-	-	1 hour
PNLIF/Com-ID Error (PNLIF to SRM)	-	-	-	-	1 hour
PNLIF/Link Up Status (PNLIF to SRM)	-	-	-	-	1 hour
PNLIF/PNLIF_Board_ID	-	-	-	-	24 hours
PNLIF/PNLIF_Board_Work	-	-	-	-	24 hours
PNLIF/Board Revision	-	-	-	-	24 hours
PNLIF/FPGA1 Revision	-	-	-	-	24 hours
PNLIF/FPGA2 Revision	-	-	-	-	24 hours
PNLIF/DSW2-1	-	-	-	-	24 hours
PNLIF/DSW2-2	-	-	-	-	24 hours
PNLIF/DSW2-3	-	-	-	-	24 hours
PNLIF/DSW2-4	-	-	-	-	24 hours
PNLIF/DSW2-5	-	-	-	-	24 hours
PNLIF/DSW2-6	-	-	-	-	24 hours
PNLIF/DSW2-7	-	-	-	-	24 hours
PNLIF/DSW2-8	-	-	-	-	24 hours
PNLIF/PNLIF Power Supply (+24V)	V	21.6 ~ 26.4	<21.6 or >26.4	-	1 hour
PNLIF/PNLIF Power Supply (+5V)	V	4.5 ~ 5.5	<4.5 or >5.5	-	1 hour
PNLIF/PNLIF Power Supply (+3.3V1)	V	2.97 ~ 3.63	<2.97 or >3.63	-	1 hour
TDRV/TDRV Board Detect Status	-	-	-	-	1 hour
TDRV/TDRV Serial Communication Error	-	-	-	-	1 hour
TDRV/Check Sum Error (SRM to TDRV)	-	-	-	-	1 hour
TDRV/Com-ID Error (SRM to TDRV)	-	-	-	-	1 hour
TDRV/Link Up Status (SRM to TDRV)	-	-	-	-	1 hour
TDRV/Check Sum Error (TDRV to SRM)	-	-	-	-	1 hour
TDRV/Com-ID Error (TDRV to SRM)	-	-	-	-	1 hour
TDRV/Link Up Status (TDRV to SRM)	-	-	-	-	1 hour
TDRV/TDRV_Board_ID	-	-	-	-	24 hours
TDRV/TDRV_Board_Work	-	-	-	-	24 hours
TDRV/Board Revision	-	-	-	-	24 hours
TDRV/FPGA1 Revision	-	-	-	-	24 hours
TDRV/DSW2-1	-	-	-	-	24 hours
TDRV/DSW2-2	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
TDRV/DSW2-3	-	-	-	-	24 hours
TDRV/DSW2-4	-	-	-	-	24 hours
TDRV/DSW2-5	-	-	-	-	24 hours
TDRV/DSW2-6	-	-	-	-	24 hours
TDRV/DSW2-7	-	-	-	-	24 hours
TDRV/DSW2-8	-	-	-	-	24 hours
TDRV/DSW3-1	-	-	-	-	24 hours
TDRV/DSW3-2	-	-	-	-	24 hours
TDRV/DSW3-3	-	-	-	-	24 hours
TDRV/DSW3-4	-	-	-	-	24 hours
TDRV/DSW3-5	-	-	-	-	24 hours
TDRV/DSW3-6	-	-	-	-	24 hours
TDRV/DSW3-7	-	-	-	-	24 hours
TDRV/DSW3-8	-	-	-	-	24 hours
TDRV/DSW4-1	-	-	-	-	24 hours
TDRV/DSW4-2	-	-	-	-	24 hours
TDRV/DSW4-3	-	-	-	-	24 hours
TDRV/DSW4-4	-	-	-	-	24 hours
TDRV/DSW4-5	-	-	-	-	24 hours
TDRV/DSW4-6	-	-	-	-	24 hours
TDRV/DSW4-7	-	-	-	-	24 hours
TDRV/DSW4-8	-	-	-	-	24 hours
TDRV/TDRV Power Supply (+24V)	V	21.6 ~ 26.4	<21.6 or >26.4	-	1 hour
TDRV/TDRV Power Supply (+5V)	V	4.5 ~ 5.5	<4.5 or >5.5	-	1 hour
TDRV/TDRV Power Supply (+3.3)	V	2.97 ~ 3.63	<2.97 or >3.63	-	1 hour
TDRV/Longitudinal Motor Voltage	V	-	-	-	1 hour
TDRV/Longitudinal Encoder Voltage	V	-	-	-	1 hour
GCSNAP/GC Snap SW CH1	-	-	-	-	1 hour
GCSNAP/GC Snap SW CH2	-	-	-	-	1 hour
GCSNAP/GC Snap SW CH3	-	-	-	-	1 hour
GCSNAP/GC Snap SW CH4	-	-	-	-	1 hour
ExtTmpHum /Shield Room Temperature	degC	5~32	<5 or >32	-	1 hour
ExtTmpHum /Shield Room Humidity	%	0~80	>80	-	1 hour
PSC/EVC LED Status (S0)	-	OFF	-	-	24 hours
PSC/EVC LED Status (S1)	-	OFF	-	-	24 hours
PSC/EVC LED Status (S2)	-	OFF	-	-	24 hours
PSC/EVC LED Status (S3)	-	OFF	-	-	24 hours
PSC/EVC LED Status (POWER ON)	-	ON	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
PSC/EVC LED Status (CONF DONE)	-	ON	-	-	24 hours
PSC/EVC LED Status (SCAN START)	-	-	-	-	24 hours
PSC/EVC LED Status (ERROR)	-	OFF	-	-	24 hours
PSC/EVC DSW2 Status	-	-	-	-	24 hours
PSC/EVC RST FPGA Revision	-	-	-	-	24 hours
PSC/EVC Board Revision	-	-	-	-	24 hours
PSC/EVC_Board_Work	-	-	-	-	24 hours
PSC/EVC_Rack_ID	-	-	-	-	24 hours
PSC/EVC_Board_ID	-	-	-	-	24 hours
PSC/EVC RF FPGA Temperature	degC	0~70	<0 or >70	-	1 hour
PSC/EVC RF FPGA Revision	-	-	-	-	24 hours
PSC/EVC GC FPGA Temperature	degC	0~70	<0 or >70	-	1 hour
PSC/EVC GC FPGA Revision	-	-	-	-	24 hours
PSC/EVC_Boot_ROM_Version	-	-	-	-	24 hours
PSC/PSC_Firmware_Version	-	-	-	-	24 hours
PSC/EVC_DSW3_Status	-	-	-	-	24 hours
DGC/DGC_Rack_ID	-	-	-	-	24 hours
DGC/DGC_Board_Work	-	-	-	-	24 hours
DGC/DGC Board Revision	-	-	-	-	24 hours
DGC/DGC PLD Revision	-	-	-	-	24 hours
DGC/DGC FPGA Revision	-	-	-	-	24 hours
DGC/DGC B0 Output Enable	-	-	-	-	24 hours
DGC/DGC GCZ Output Enable	-	-	-	-	24 hours
DGC/DGC GCY Output Enable	-	-	-	-	24 hours
DGC/DGC GCX Output Enable	-	-	-	-	24 hours
DGC/DGC RMFC Enable	-	-	-	-	24 hours
DGC/DGC VFC Output Enable	-	-	-	-	24 hours
DGC/DGC ECC Output Enable	-	-	-	-	24 hours
DGC/DGC DSW1-1	-	-	-	-	24 hours
DGC/DGC DSW1-2	-	-	-	-	24 hours
DGC/DGC DSW1-3	-	-	-	-	24 hours
DGC/DGC DSW1-4	-	-	-	-	24 hours
DGC/DGC DSW2-1	-	-	-	-	24 hours
DGC/DGC DSW2-2	-	-	-	-	24 hours
DGC/DGC DSW2-3	-	-	-	-	24 hours
DGC/DGC DSW2-4	-	-	-	-	24 hours
DGC/DGC CNT FPGA Temperature	degC	0~70	<0 or >70	-	1 hour
DGC/DGC LED (Software)	-	OFF	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
DGC/DGC LED (System Clock)	-	-	-	-	24 hours
DGC/DGC LED (DGC Error)	-	OFF	-	-	24 hours
DGC/DGC LED (Mode)	-	OFF	-	-	24 hours
DGC/DGC LED (GCZ Enable)	-	ON	-	-	24 hours
DGC/DGC LED (GCY Enable)	-	ON	-	-	24 hours
DGC/DGC LED (GCX Enable)	-	ON	-	-	24 hours
DGC/DGC LED (RFGC Active)	-	-	-	-	24 hours
DGC/Power Supply (3.3VD)	-	3.3	-	-	24 hours
DGC/Power Supply (1.1V)	-	1.1	-	-	24 hours
DGC/Power Supply (2.5V)	-	2.5	-	-	24 hours
DGC/Power Supply (3.3V)	-	3.3	-	-	24 hours
DGC/Power Supply (12VD)	-	12	-	-	24 hours
DGC/Command Overflow (GZ)	-	Normal	Warning	-	24 hours
DGC/Command Overflow (GY)	-	Normal	Warning	-	24 hours
DGC/Command Overflow (GX)	-	Normal	Warning	-	24 hours
DGC/Improper Command Warning (GZ)	-	Normal	Warning	-	24 hours
DGC/Improper Command Warning (GY)	-	Normal	Warning	-	24 hours
DGC/Improper Command Warning (GX)	-	Normal	Warning	-	24 hours
DGC/Operating Clock Warning	-	Normal	Warning	-	24 hours
DGC/SHIMAMP Error	-	Normal	-	Fault	24 hours
DGC/SHIMAMP I/F Error	-	Normal	-	Fault	24 hours
DGC/Improper Command Error (GZ)	-	Normal	-	Fault	24 hours
DGC/Improper Command Error (GY)	-	Normal	-	Fault	24 hours
DGC/Improper Command Error (GX)	-	Normal	-	Fault	24 hours
DGC/GC Serial Command Link Error	-	Normal	-	Fault	24 hours
DGC/DGC_Board_ID	-	-	-	-	24 hours
DGC/SHIMAMP BRD1CH1 Enable	-	-	-	-	24 hours
DGC/SHIMAMP BRD1CH2 Enable	-	-	-	-	24 hours
DGC/SHIMAMP BRD1CH1 ErrorStatus	-	-	-	-	24 hours
DGC/SHIMAMP BRD1CH2 ErrorStatus	-	-	-	-	24 hours
DGC/SHIMAMP Interface ErrorStatus	-	-	-	-	24 hours
DGC/SHIMAMP B1CH1 Enable	-	-	-	-	24 hours
DGC/SHIMAMP B1CH1 OverCurrent	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 OverTemp	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 OpenLoad	-	-	-	Fault	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
DGC/SHIMAMP B1CH1 EXT-CLKWarning	-	-	Warning	-	24 hours
DGC/SHIMAMP B1CH1 Over DCBusVol	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 Desat Fault	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 Output Oscillation	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 PLL Fault	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 Com Timeout	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH1 Checksum Fault	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 Enable	-	-	-	-	24 hours
DGC/SHIMAMP B1CH2 OverCurrent	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 OverTemp	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 OpenLoad	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 EXT-CLKWarning	-	-	Warning	-	24 hours
DGC/SHIMAMP B1CH2 Over DCBusVol	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 Desat Fault	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 Output Oscillation	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 PLL Fault	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 Com Timeout	-	-	-	Fault	24 hours
DGC/SHIMAMP B1CH2 Checksum Fault	-	-	-	Fault	24 hours
COILCONT/COILCONT_Board_ID	-	-	-	-	24 hours
COILCONT/COILCONT_Board_Work	-	-	-	-	24 hours
COILCONT/Board Revision	-	-	-	-	24 hours
COILCONT/FPGA1 Revision	-	-	-	-	24 hours
COILCONT/FPGA2 Revision	-	-	-	-	24 hours
COILCONT/FPGA3 Revision	-	-	-	-	24 hours
COILCONT/FPGA4 Revision	-	-	-	-	24 hours
COILCONT/DIPSW Status	-	-	-	-	24 hours
COILCONT/COILCONT Module Error Status	-	Normal	-	-	1 hour
COILCONT/COILIF PCB Error Status	-	Normal	-	-	1 hour

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
COILCONT/CN1 COILDRV PCB Part Error	-	Normal	-	-	1 hour
COILCONT/CN2 COILDRV PCB Part Error	-	Normal	-	-	1 hour
COILCONT/CN3 COILDRV PCB Part Error	-	Normal	-	-	1 hour
COILCONT/CN4 COILDRV PCB Part Error	-	Normal	-	-	1 hour
COILCONT/Tx COILDRV PCB Part Error	-	Normal	-	-	1 hour
COILCONT/COILCONT - SYNTH COM	-	Normal	-	-	1 hour
COILCONT/COILCONT - LINESEL COM	-	Normal	-	-	1 hour
COILCONT/FPGA COM (MASTER - SLAVE1)	-	Normal	-	-	1 hour
COILCONT/FPGA COM (MASTER - SLAVE2)	-	Normal	-	-	1 hour
COILCONT/FPGA COM (MASTER - SLAVE3)	-	Normal	-	-	1 hour
COILCONT/CN1 COIL ID COM	-	Normal	-	-	1 hour
COILCONT/CN2 COIL ID COM	-	Normal	-	-	1 hour
COILCONT/CN3 COIL ID COM	-	Normal	-	-	1 hour
COILCONT/CN4 COIL ID COM	-	Normal	-	-	1 hour
COILCONT/Input Power Supply (+35V)	-	Normal	-	Error	1 hour
COILCONT/Input Power Supply (-23V)	-	Normal	-	Error	1 hour
COILCONT/Input Power Supply (+26V)	-	Normal	-	Error	1 hour
COILCONT/Input Power Supply (-14.5V)	-	Normal	-	Error	1 hour
SYNTHE/DSF FPGA Temperature	-	0~70	<0 or >70	-	1 hour
SYNTHE/DSF LED Status (Reserve1)	-	OFF	-	-	24 hours
SYNTHE/DSF LED Status (DEC ERR)	-	OFF	-	-	24 hours
SYNTHE/DSF LED Status (Reserve2)	-	OFF	-	-	24 hours
SYNTHE/DSF LED Status (DBG MODE)	-	-	-	-	24 hours
SYNTHE/DSF LED Status (PLL LOCK)	-	ON	-	-	24 hours
SYNTHE/DSF LED Status (RFPA ON)	-	ON	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
SYNTHE/DSF LED Status (RFPA BIAS)	-	-	-	-	24 hours
SYNTHE/DSF LED Status (RESET)	-	-	-	-	24 hours
SYNTHE/DSF CLKRST FPGA	-	-	-	-	24 hours
SYNTHE/DSF FPGA Revision	-	-	-	-	24 hours
SYNTHE/DSF Frequency ROM Revision	-	-	-	-	24 hours
SYNTHE/DSF Frequency ROM Offset	-	-	-	-	24 hours
SYNTHE/DSF Frequency ROM Hardware	-	-	-	-	24 hours
SYNTHE/DSF DDS1 Channel A Initial Phase Value	-	-	-	-	24 hours
SYNTHE/DSF DDS1 Channel B Initial Phase Value	-	-	-	-	24 hours
SYNTHE/DSF DDS2 Channel A Initial Phase Value	-	-	-	-	24 hours
SYNTHE/DSF DDS2 Channel B Initial Phase Value	-	-	-	-	24 hours
SYNTHE/DSF_Rack_ID	-	-	-	-	24 hours
SYNTHE/DSF_Board_ID	-	-	-	-	24 hours
SYNTHE/DSF Board Revision	-	-	-	-	24 hours
SYNTHE/DSF_Board_Work	-	-	-	-	24 hours
EXIF/EXIF_Rack_ID	-	-	-	-	24 hours
EXIF/EXIF_Board_ID	-	-	-	-	24 hours
EXIF/EXIF_Board_Work	-	-	-	-	24 hours
EXIF/EXIF Board Revision	-	-	-	-	24 hours
EXIF/EXIF_FPGA_Revision	-	-	-	-	24 hours
EXIF/PMM Select ID	-	-	-	-	24 hours
EXIF/EXIF LED (S0)	-	-	-	-	24 hours
EXIF/EXIF LED (S1)	-	-	-	-	24 hours
EXIF/EXIF LED (S2)	-	-	-	-	24 hours
EXIF/EXIF LED (S3)	-	-	-	-	24 hours
EXIF/EXIF LED (POWER)	-	-	-	-	24 hours
EXIF/EXIF LED (DAT Error)	-	-	-	-	24 hours
EXIF/EXIF LED (CPU IRQ)	-	-	-	-	24 hours
EXIF/EXIF LED (PMM Trigger)	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-1	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-2	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
EXIF/EXIF DIPSW2-3	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-4	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-5	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-6	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-7	-	-	-	-	24 hours
EXIF/EXIF DIPSW2-8	-	-	-	-	24 hours
EXIF/EXIF FPGA Temperature	degC	0~70	<0 or >70	-	1 hour
EXIF/EXIF RST FPGA Revision	-	-	-	-	24 hours
EXIF/EXIF_Boot_ROM_Version	-	-	-	-	24 hours
EXIF/EXIF_Firmware_Version	-	-	-	-	24 hours
ADIF/ADIF_Board_ID	-	-	-	-	24 hours
ADIF/ADIF_Board_Work	-	-	-	-	24 hours
ADIF/ADIF Board Revision	-	-	-	-	24 hours
ADIF/ADIF RST FPGA Revision	-	-	-	-	24 hours
ADIF/ADIF FPGA Revision	-	-	-	-	24 hours
ADIF/ADIF LED (Reserved)	-	-	-	-	24 hours
DIF/ADIF LED (Power ON(12V))	-	-	-	-	24 hours
ADIF/ADIF LED (Debug AD Mode)	-	-	-	-	24 hours
ADIF/ADIF LED (Debug Clock Mode)	-	-	-	-	24 hours
ADIF/ADIF LED (FPGA PLL Locked)	-	-	-	-	24 hours
ADIF/ADIF LED (Error)	-	-	-	-	24 hours
ADIF/ADIF LED (ADC Start)	-	-	-	-	24 hours
ADIF/ADIF LED (Receiving Mode)	-	-	-	-	24 hours
ADIF/ADIF LED (CONF DONE)	-	-	-	-	24 hours
ADIF/ADIF LED (CONF ERR)	-	-	-	-	24 hours
ADIF/ADIF LED (PWR OK)	-	-	-	-	24 hours
ADIF/ADIF LED (MST RST)	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-1	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-2	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-3	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-4	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-5	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-6	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-7	-	-	-	-	24 hours
ADIF/ADIF DIPSW2-8	-	-	-	-	24 hours
ADIF/ADIF FPGA Temperature Reg	degC	0~70	<0 or >70	-	1 hour
ADIF/ADIF_Rack_ID	-	-	-	-	24 hours
DSR1/DSR1_Board_ID	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
DSR1/DSR1 Board Work	-	-	-	-	24 hours
DSR1/DSR1 Board Revision	-	-	-	-	24 hours
DSR1/DSR1 RST FPGA Revision	-	-	-	-	24 hours
DSR1/DSR1 FPGA Revision	-	-	-	-	24 hours
DSR1/DSR1 LED (AD TRIG)	-	-	-	-	24 hours
DSR1/DSR1 LED (ADC LINK)	-	-	-	-	24 hours
DSR1/DSR1 LED (POWER ON)	-	-	-	-	24 hours
DSR1/DSR1 LED (CONF DONE)	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-1	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-2	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-3	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-4	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-5	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-6	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-7	-	-	-	-	24 hours
DSR1/DSR1 DIPSW2-8	-	-	-	-	24 hours
DSR1/RECEIVE1 FPGA Temperature Reg	degC	0~70	<0 or >70	-	1 hour
DSR1/DSR1 Rack ID	-	-	-	-	24 hours
DSR1/DSR1 Voltage (3.3V) Mon1 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (2.5V) Mon2 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA1.1V) Mon3 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA1.15V) Mon4 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA1.8V) Mon5 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA2.5V) Mon6 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA3.3V) Mon7 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (FPGA1.5V) Mon8 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (None) Mon9 Reg	-	-	-	-	24 hours
DSR1/DSR1 Voltage (5V) Mon10 Reg	-	-	-	-	24 hours
DSR2/DSR2_Board_ID	-	-	-	-	24 hours
DSR2/DSR2_Board_Work	-	-	-	-	24 hours
DSR2/DSR2 Board Revision	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
DSR2/DSR2 RST FPGA Revision	-	-	-	-	24 hours
DSR2/DSR2 FPGA Revision	-	-	-	-	24 hours
DSR2/DSR2 LED (AD TRIG)	-	-	-	-	24 hours
DSR2/DSR2 LED (ADC LINK)	-	-	-	-	24 hours
DSR2/DSR2 LED (POWER ON)	-	-	-	-	24 hours
DSR2/DSR2 LED (CONF DONE)	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-1	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-2	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-3	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-4	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-5	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-6	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-7	-	-	-	-	24 hours
DSR2/DSR2 DIPSW2-8	-	-	-	-	24 hours
DSR2/RECEIVE2 FPGA Temperature Reg	degC	0~70	<0 or >70	-	1 hour
DSR2/DSR2_Rack_ID	-	-	-	-	24 hours
DSR2/DSR2 Voltage (3.3V) Mon1 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (2.5V) Mon2 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA1.1V) Mon3 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA1.15V) Mon4 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA1.8V) Mon5 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA2.5V) Mon6 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA3.3V) Mon7 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (FPGA1.5V) Mon8 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (None) Mon9 Reg	-	-	-	-	24 hours
DSR2/DSR2 Voltage (5V) Mon10 Reg	-	-	-	-	24 hours
Common					
Monitor Mode	-	OFF	-	-	1 hour
BE PC/CPU Model	-	-	-	-	24 hours
BE PC/Actual Physical Memory	byte	-	< 17179869184 (MR-IRCP-M1)	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
BE PC/Total Physical Memory	byte	-	≠8446652416(MR-RFIP-1N) ≠16996012032(MR-IRCP-M1)	-	24 hours
BE PC/Available Physical Memory	KB	131072~	-	<131072	1 hour
BE PC/Mother Board	-	-	-	-	24 hours
BE PC/Monitor Manufacturer	-	-	-	-	24 hours
BE PC/DVD-Multi Drive	-	-	-	-	24 hours
BE PC/Disk Drive	-	-	-	-	24 hours
BE PC/HDD Size C:	byte	-	-	-	24 hours
BE PC/HDD Size Q:	byte	-	-	-	24 hours
MR-IRCP-M1					
BE PC/HDD Size P:	byte	-	-	-	24 hours
BE PC/HDD Size R:	byte	-	-	-	24 hours
Common					
BE PC/HDD Free Space C:	byte	1073741824~	-	<1073741824	1 hour
BE PC/HDD Free Space Q:	byte	-	-	-	1 hour
MR-IRCP-M1					
BE PC/HDD Free Space P:	byte	-	-	-	1 hour
BE PC/HDD Free Space R:	byte	-	-	-	1 hour
Common					
BE PC/HDD Life Time	hour	-	-	-	1 hour
BE PC/HDD Temperature	degC	10~55	<10 or (55< & <=60)	>60	1 hour
BE PC/BIOS Version	-	-	-	-	24 hours
BE PC/App Mem Usage	byte	-	-	-	1 hour
BE PC/IP Mem Usage	byte	-	-	-	1 hour
BE PC/Database Size	byte	~3650722200	3650722200< & <=4037269258	>4037269258	1 hour
BE PC/DB Transaction Size	byte	-	-	-	1 hour
BE PC/Start Stop Count	-	0~50000	>50000	-	1 hour
BE PC/Current Pending Sector Count	-	-	-	-	1 hour
BE PC/Offline Uncorrectable	-	0	>0	-	1 hour
BE PC/Reallocated Sector Ct	-	-	-	-	1 hour
BE PC/UltraDMA CRC Error Count	-	-	-	-	1 hour
BE PC/Mother Voltage	V	2.6~3.2	<2.5	-	1 hour
FE PC/CPU Model	-	-	-	-	24 hours
FE PC/Actual Physical Memory	byte	-	< 34359738368 (MR-IRCP-M1)	-	24 hours
FE PC/Total Physical Memory	byte	-	≠6123008000(MR-RFIP-1N) ≠26380460032(MR-IRCP-M1)	-	24 hours
FE PC/Available Physical Memory	KB	131072~	<131072	-	1hour
FE PC/Mother Board	-	-	-	-	24 hours

Items	Unit	Normal range	Alert conditions	Error conditions	Interval
FE PC/DVD-ROM Drive	-	-	-	-	24 hours
FE PC/Disk Drive	-	-	-	-	24 hours
FE PC/HDD Size C:	byte	-	-	-	24 hours
FE PC/HDD Size D:	byte	-	-	-	24 hours
FE PC/HDD Free Space C:	byte	536870912~	<536870912	-	1 hour
FE PC/HDD Free Space D:	byte	536870912~	<536870912	-	1 hour
FE PC/HDD Life Time	hour	-	-	-	1 hour
FE PC/HDD Temperature	degC	10~55	<10 or (55< & <=60)	>60	1 hour
FE PC/BIOS Version	-	-	-	-	24 hours
FE PC/Start/Stop Count	-	0~50000	>50000	-	1 hour
FE PC/Reallocated Sector Ct	-	-	-	-	1 hour
FE PC/Current Pending Sector Count	-	-	-	-	1 hour
FE PC/Offline Uncorrectable	-	0	>0	-	1 hour
FE PC/UltraDMA CRC Error Count	-	-	-	-	1 hour
FE PC/Mother Voltage	V	2.6~3.2	<2.5	-	1 hour
Work History/Change Cold Head Date	-	-	-	-	-
Work History/Change HDD Date BE	-	-	-	-	-
Work History/Change HDD Date FE	-	-	-	-	-
Work History/Change Mother Battery Date BE	-	-	-	-	-
Work History/Change Mother Battery Date FE	-	-	-	-	-

NOTE Total Physical Memory is assigned value.

In normal state, the value should be zero. Other value indicates abnormal state.

7.2.11 Monitoring event log (system) errors.....

ECHELON Smart monitors event log of Windows. If it detects errors, it notifies of that by e-mail.

If it detects errors, execute “Collect Log” from [SERVICE] - [ACTIONS] - [Run command], and refer to event log. (Refer to Chapter7.4.2.2.)

Table7-8 Sources of Error Event

#	Sources of Error Event	The main cause of errors
1	adpu160m	The host adapter has shown abnormalities.
2	Alerter	
3	AsyncMac	
4	atapi	The drive has shown abnormalities.
5	Atmarpc	
6	AtmElan	
7	ATMEPVCP	
8	AtmLane	
9	Atmuni	
10	Browser	The browser has shown abnormalities.
11	Cdm	
12	Dhcp	In the case of share violations, DHCP server stops leasing addresses.
13	dmbboot	Disks can not be read.
14	dmio	Disks can not be read or written.
15	Dnsapi	DNS can not be registered.
16	Dnscache	
17	E1000	LAN cards or hubs have shown abnormalities.
18	ETG2P51	
19	fdisk	Bad sectors are detected in the disk.
20	Http	Permission of the directory has been changed.
21	i8042prt	The mouse has shown abnormalities.
22	IAS	Authentication of the Internet fails.
23	IGMPv2	
24	IISCTLS	
25	IISLOG	Access to the folder is restricted.
26	IISMAP	
27	intelppm	
28	IPBOOTP	
29	IPMGM	

#	Sources of Error Event	The main cause of errors
30	IPNATHLP	The host exists in the network where IP addresses are duplicated in address NAT DHCP pool.
31	IPRIP2	Internet Control Message Protocol (ICMP) can not be reached.
32	IPRouterManager	No network adapter is detected.
33	IPSec	The policy store has damaged files.
34	IPXCP	
35	IPXRIP	
36	IPXRouterManager	00000000 is set as the number of the internal network of NWLink IPX/SPX protocol installed in the computer.
37	IPXSAP	
38	isapnp	
39	kbdclass	Drag of mouse operation failed.
40	LDM	Disks can not be read or written.
41	LDMS	Disks can not be read or written.
42	LmHosts	
43	mouclass	No mouse is recognized.
44	mouhid	
45	MRxDAV	
46	MrxSmb	Computer names are duplicated in the network.
47	MSFTPSVC	IPSec is not removed from MetaBase correctly after IP restriction was changed.
48	NdisWan	
49	NetBIOS	
50	NetBT	Initializing NetBT fails.
51	Netlogon	Domain authentication failed.
52	Nla	
53	NWCWorkstation	
54	NWRDR	The network has shown abnormalities.
55	OSPF	
56	OSPFMib	
57	Parport	
58	PolicyAgent	IPSec connection has shown abnormalities.
59	PptpMiniport	
60	PSched	The network adapter has shown abnormalities.
61	RasAuto	
62	RasMan	
63	Rawwan	
64	Rdbss	

#	Sources of Error Event	The main cause of errors
65	redbook	The driver has shown abnormalities.
66	RemoteAccess	IP address setting has shown abnormalities.
67	Removable Storage Service	RSM can not automatically configure robotic libraries.
68	Schannel	
69	Serial	No registry key exists.
70	Server	There is the share folder that does not exist.
71	Srv	Files in the network can not be read.
72	StillImage	
73	Tcpip	Network connection has shown abnormalities.
74	Tcpip6	
75	TCPMon	
76	tunmp	
77	Update	
78	UPS	
79	vncdrv	
80	VolSnap	The volume of difference area can not be used after boot.
81	W32Time	Time synchronization fails.
82	W3SVC	Files can not be accessed.
83	WAM	Mismatched dynamic link libraries exist.
84	Windows Script Host	Abnormalities occur in the application of Windows.
85	Workstation	WINS database has static entry.
86	WZCSVC	

7.2.12 Monitor Service.....

If a customer chooses “Monitor Service”, the system is monitored although [GANTRY] button is being turned off. If an abnormality is detected, it is reported by e-mail.




Figure7-31 Mail notification

Table7-9 Monitoring items at the time of Monitor Service

#	Items	Error condition
1	Magnet System(Alarm)	When the "Magnet System(Alarm)" is Fault.
2	Magnet System(Warning)	When the "Magnet System(Warning)" is Warning.
3	Magnet System (Quench Alarm)	When the "Magnet System (Quench Alarm)" is Fault.
4	Compressor System(Alarm)	When the "Compressor System(Alarm)" is Warning.
5	Magnet Critical Error	When the "Magnet Critical Error" is Fault.
6	Compressor Flow Alarm	When the "Compressor Flow Alarm" is Fault.
7	Primary Supply Temp Alarm	When the "Primary Supply Temp Alarm" is Fault.

7.2.13 Tool & Manual.....


By Tool & Manual, Sentinel Customer Support manual and related tools can be acquired.

The following is the procedure.


- 1 Select [GROUPS] - [Tool & manual].




- 2 Click [SERVICE] - [FILE MANAGEMENT]-[Manage system files].




3 System Files screen is displayed.



4 Select the desirable file to download and click "Download". (Multiple selection is accepted.)




- 5 The following window is displayed, so click a link.




- 6 The dialog of file download is displayed, so click "Save".



7 Specify saving place, the file name and click [Save].



8 [Download complete] dialog is displayed, so click [Close] button.



7.3 Manage alerts

[Manage alerts] indicates Alerts occurring in device and enables [Acknowledge Alerts] and [Close Alerts]. Confirm Alert detail through Web and perform [Acknowledge Alerts] and [Closed Alerts] if receiving Alert mail.




Figure7-32 Flowchart from Alert [Acknowledge Alerts] to [Closed Alerts].




Figure7-33 Alerts state transition diagram

Table7-10 Alert status

Definitions of alert states		
Active	Pending	This is the alert request's initial state when it is created. The alert request is "pending" processing.
	Submitted	The processing of the alert request has begun; it has been submitted to the alert processor.
	Sent	The processing of the alert request is completed; it has been successfully sent to the alert processor.
	Acknowledged	A user has reviewed the alert request and acknowledged that it exists.
Completed	Closed	The alarm condition has been cleared or otherwise handled, and no more user action is required. The alert request is closed.
	Duplicate	The alert request pertains to the same alarm and alert as another alert request for the same device that is currently being processed. This alert request, therefore, is a duplicate and can be ignored.
	Suppressed	The alert escalator has exhausted all escalation levels for the alert request, and processing has reached the end. Therefore, all automatic processing is now suppressed.

If you select any group or device and choose [SERVICE] - [ACTIONS] - [Manage alerts], alert managing screen is displayed.




Figure7-34 Manage alerts

Alert description is displayed as follows. Alert is acknowledged and closed here.





Figure7-35 Alerts description screen

- (1) Alert
Alert name is indicated.
- (2) Status
Alert status name is indicated.
- (3) Device
Alert occurrence device name is indicated.
NOTE Nothing is indicated if device is selected.
- (4) Date
Alert occurrence date is indicated.
- (5) Switching screen to display.
By clicking this tab, screen is switched to filter screen.
- (6) Alert acknowledgement.
By clicking this button, selected alert is acknowledged.
- (7) Help
By clicking , online help is displayed.
- (8) Close
By clicking [Close] button, screen closes.

If you click **FILTER** tab, screen is switched to filter screen and it becomes possible to filter alert description with status, date to search.





Figure7-36 Alert filter screen

- (1) Status
Alert status is selected here.
- (2) Start date
Date range specification start date is filled in.
- (3) End date
Date range specification end date is filled in.
- (4) Search
If clicking **Search** , search is performed according to specified filter.
- (5) Display screen switch
If clicking this tab, screen is switched to Alert screen.
- (6) Help
By clicking **?** , online help is displayed.
- (7) Close
by clicking [Close] button, this screen closes.

The following is search procedure.

- 1 Select Alert status to display. "All" is selected by default.




- 2 If specifying date, fill in start date and end date. Entry order is "Month/Day/Year Time:Minutes"

NOTE It is not possible to filter date unless month/day/year is filled in.


Start date / /  :

End date / /  :

Click  and calendar is displayed.




Date can be specified from calendar.






- (1) Month selection
Displayed calendar month is selected.
- (2) Year selection
Displayed calendar year is selected.
- (3) Cancel
Specifying date from calendar is canceled.
- (4) Date selection
Specification date is clicked.

Calendar indicates execution date month. If you click ◀, the previous calendar is displayed instead. If you click ▶, coming month calendar is displayed instead.


Calendar can be specified from year and month.



If clicking date to specify from calendar, date can be filled in.





- 3 After filtering, click **Search** and specified Alert description is indicated.




The Alerts window shows a search interface. The 'Status' dropdown is set to 'Suppressed'. The 'Start date' is '1 / 1 / 2006' and the 'End date' is '1 / 31 / 2006'. A red box highlights the 'Search' button.

Field	Value
Status	Suppressed
Start date	1 / 1 / 2006
End date	1 / 31 / 2006




The Alerts window displays a table with one search result. The table has columns for Alert, Status, Device, and Date.

Alert	Status	Device	Date
 Error_Alert	Suppressed	MRI_ECHELON_Smart...	2006/01/23 14:38:19


The following is Alert acknowledging procedure.


- 1 To acknowledge alert, choose alert to acknowledge.




To select Alerts, by using [Ctrl] key and [Shift] key, multiple selections are possible, too.



- 2 Click after selecting alert to acknowledge. If clicking , alert acknowledging screen is displayed.





- (1) Alert
Selected alert name is indicated.
- (2) Alarm
Selected alert alarm name is indicated.
- (3) Status
Selected alert status type is indicated.
- (4) Date
Alert occurrence date is indicated.
- (5) Assign to
Approver is selected.
NOTE Only administrator user can use this function.
- (6) Notes
Alert status change reason is entered.
- (7) Acknowledge Alerts
If you click , alert is acknowledged.
- (8) Close Alerts
If you click , alert is closed.
- (9) Cancel
If you click , alert acknowledgement is canceled and this screen is closed.
- (10) Help
If you click  , online help is displayed.




- 3 Click Acknowledge Alerts to acknowledge alert.



- 4 Alert status is "Acknowledged" if it is acknowledged.

The following is alert closing procedure.

- 1 Click Close Alerts to close alert.



- 2 If alert is closed, 🚩 is not marked any more and its status becomes "Closed".

7.4 Run command

7.4.1 Basic operation

Device command can be run. Select device and then click [SERVICE] - [ACTIONS] - [Run command].





Figure7-37 Run command




Figure7-38 COMMANDS

- (1) Command name
Running command is selected.
- (2) Screen switch
[PARAMETERS] screen and [COMMAND RESULTS] screen are switched.
- (3) Run
If you click , selected command is run.
- (4) Close
If you click , this screen closes.
- (5) Help
If you click

If you open [COMMAND RESULTS] tab, next screen is displayed.




Figure7-39 COMMAND RESULTS

- (1) Command
Command name is indicated.
- (2) Description
Command description is indicated.
- (3) Start Date
Command running date is indicated.
- (4) End Date
Command ending date is indicated.
- (5) Status
Executed command status is indicated.
- (6) By
Name of the user who ran command is indicated.
- (7) INPUT PARAMETERS
If command running result is selected, input parameters are indicated.
NOTE Nothing is indicated at the time of no parameter.
- (8) OUTPUT PARAMETERS
If command running result is selected, output parameter is indicated.

7.4.2 ECHELON Smart command.....

ECHELON Smart command has the following functions.

- Screen Capture See 7.4.2.1.
- Collect Log See 7.4.2.2.
- Service Backup See 7.4.2.3.
- Service Restore See 7.4.2.4.
- Disk Cleanup See 7.4.2.5.
- Shutdown See 7.4.2.6.
- Collect Manage Support Log See 7.4.2.7.
- Get Prescan Data See 7.4.2.8.
- Password Reset See 7.4.2.9.
- Collect Hardware Status Log See 7.4.2.10.

7.4.2.1 Screen Capture


Screen Capture is possible in order to confirm screen status at the time of failure occurrence.

NOTE Device screen automatically switches to personal information protection mode at the time of Screen Capture, so no individual information indicated on screen is captured.
See chapter 7.7.4.


Acquire Screen Capture according to the following procedure.

(1) Screen Capture execution

- 1 Select Screen Capture target device and click [Run command].



2 Select [Screen Capture] and click [Run] button.




3 Open [COMMAND RESULTS] tab and executed command "Status" is "Started".

Command	Description	Start Date	End Date	Status	By
CaptScrn	Screen Capture	2006/09/1...		Started	hirota

4 Close the command screen once, take 1 or 2 minutes and then click [Refresh] button at the upper right of the screen.

- 5 Click [Run command] again and open [COMMAND RESULTS] tab.



- 6 Status indicates "Completed" if process finishes. Select the command to have run.

This is a close-up view of the 'COMMAND RESULTS' table from the previous screenshot. Red boxes highlight the 'CaptScr' command in the first row and the 'Completed' status in the same row.

Command	Description	Start Date	End Date	Status	By
CaptScr	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota

- 7 [OUTPUT PARAMETERS] value indicates “ok” if process is normally completed.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
CaptScrn	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota
		2006/09/1...	2006/09/1...	Completed	hirota


INPUT PARAMETERS		OUTPUT PARAMETERS	
Parameter	Value	Parameter	Value
		RtnCaptScrn	ok

- 8 Screen Capture file is automatically uploaded to server in a few minutes.


(2) How to download and display from server

If downloading Screen Capture from server to user's PC, operate as follows.


- 1 Select the device where Screen Capture was carried out and click [SERVICE] - [FILE MANAGEMENT]-[Manage system files].



- 2 The screen of System Files is displayed, so select “Screen_{the name of the system}_{“Year_month_day_hour_minute_second”}.png” and click “Download”.




- 3 The following window is displayed, so click a link.




- 4 The dialog of file download is displayed, so click “Save”.



- 5 Specify saving place, the file name and click [Save].



- 6 [Download complete] dialog is displayed, so click [Close] button.




7.4.2.2 Collect Log

To investigate device at the time of system down and failure, system log can be acquired.


Acquire system log according to the following procedure.

(1) Collect Log execution

- 1 Select the device for Collect Log and click [Run command].



- 2 Select [Collect Log] and click [Run].



- 3 Open [COMMAND RESULTS] tab and executed command "Status" is "Started".

Command	Description	Start Date	End Date	Status	By
GetLog	Collect Log	2006/09/1...		Started	hirota

- 4 Close COMMANDS screen and click [Refresh] upper-right corner of the screen.

NOTE Take five minutes before updating.

- 5 Click [Run command] again and open [COMMAND RESULTS] tab.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
GetLog	Collect Log	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Timed out	hirota
CaptScrn	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota
		2006/09/1...	2006/09/1...	Completed	hirota

INPUT PARAMETERS		OUTPUT PARAMETERS	
Parameter	Value	Parameter	Value
		RtnGetLog	ok

- 6 “Completed” character is displayed in “Status” when process is completed. Select run commands.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
GetLog	Collect Log	2006/09/1...	2006/09/1...	Completed	hirota

- 7 “ok” is displayed in the value of [OUTPUT PARAMETER] when process is normally completed.


OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	ok

- 8 The log file of "diag_dump_info_{System name}_{Serial number}_{"MonthDay"}_{"HourMinuteSecond"}.zip" is automatically uploaded into the server in a few minutes. "BELogFile.zip", "ScanLogFile.zip" and "ReconLogFile.zip" are uploaded at the time of system reboot.


(2) Way of downloading from server and displaying

Download Log to pc you use from server at the time of log file reference. Operate according to the following procedure.


- 1 Select the device for Collect Log and click [SERVICE] - [FILE MANAGEMENT] - [Manage system files].




- The screen of System Files is displayed, so select “diag_dump_info_{System name}_{Serial number}_{"MonthDay"}_{"HourMinuteSecond"}.zip” and click “Download”.




- The following window is displayed, so click a link.




- 4 The dialog of file download is displayed, so click “Save”.



- 5 Specify saving place, the file name and click [Save].



- 6 [Download complete] dialog is displayed, so click [Close] button.




- 7 Execute the same operation to "BELogFile.zip", "ScanLogFile.zip" and "ReconLogFile.zip", too.

7.4.2.3 Service Backup


Device service back up can be acquired by remote in case of system down.

Get service backup file according to the following procedure.

- 1 Select service backup acquisition target device and click [Run command].



- 2 Select [Service Backup] and click [Run].



- 3 Open [COMMAND RESULTS] tab and executed command “Status” is “Started”.

Command	Description	Start Date	End Date	Status	By
GetSB	Service Backup	2006/09/1...		Started	hirota

- 4 Close COMMANDS screen once and click [Refresh] upper-right corner of the screen.

NOTE Take five minutes before updating.

- 5 Click [Run command] again and open [COMMAND RESULTS] tab.

Command	Description	Start Date	End Date	Status	By
GetSB	Service Backup	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Timed out	hirota
CaptScrn	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota
		2006/09/1...	2006/09/1...	Completed	hirota

INPUT PARAMETERS		OUTPUT PARAMETERS	
Parameter	Value	Parameter	Value
		RtnGetSB	ok

- 6 “Completed” character is displayed in “Status” when process is completed. Select run commands.

Command	Description	Start Date	End Date	Status	By
GetSB	Service Backup	2006/09/1...	2006/09/1...	Completed	hirota

- 7 “ok” is displayed in the value of [OUTPUT PARAMETERS] when process is normally completed.

OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	ok

- 8 The service backup file is automatically uploaded on the server after a few minutes.

7.4.2.4 Service Restore


The restoration to the device needs to get the service backup. At this time, server saving service backup file needs to be transferred to the device in advance. Since this function uses the [File transfer] function, execute the procedure of 7.6.1 Device file information update in advance. Perform transfer action (download) to the device according to the following procedure.


(1) Transfer a backup file from the server to the device


- 1 Before executing Service Restore command, ServiceBackup data to send MRI system have to be prepared.

※The ServiceBackup data to send MRI system is in user's PC.

Click [Manage software packages] of [CONFIGURATION]—[SYSTEM CONFIGURATION].



2 Click  button on Software Packages dialog.



- 3 Input Name and Description on Software Package dialog. The Name and Description will be displayed on Software Packages dialog.

Input Install Command as follows.

C:\Program Files\SENTINEL\bin\SPMWrapper.vbs

Software Package

Name: ServiceBackup
Description: ServiceBackup Data

Prerequisite Configuration

- System Software Information
- System Software Information (Echelon)
- HMC
- HMSA
- System Software Information (OASIS)
- System Software Information (XR)

Resulting configuration: -- Select --

Install Type

Automatic Manual

Locked

Package File Name: +
ServiceBackup.zip

File	Description	Last Modified
ServiceBackup.zip		Jan 17, 20...
ServiceBackup.zip		Jan 18, 20...
ServiceBackup.zip		Jan 30, 20...

COMMANDS

Pre-download command:

Post-download command:

Pre-install command:

Install Command: C:\Program Files\SENTINEL\bin\SPMWrapper.vbs

Recover Command:

Post-install command:


Pre-activate command:

DEVICE TYPES

Name	Manufacturer	Product	Model	Version
<input type="checkbox"/> AIRIS Elite	Hitachi Medi...	MRI	AIRIS Elite	V5.0
<input type="checkbox"/> AIRIS Elite...	Hitachi Medi...	MRI	AIRIS Elite...	V1.0
<input type="checkbox"/> AIRIS Elite...	Hitachi Medi...	MRI	AIRIS Elite...	V1.0
<input type="checkbox"/> AIRIS Light	Hitachi Medi...	MRI	AIRIS Light	V1.0
<input type="checkbox"/> AIRIS Light...	Hitachi Medi...	MRI	AIRIS Light...	V1.0
<input type="checkbox"/> AIRIS Soleil	Hitachi Medi...	MRI	AIRIS Soleil	V1.0
<input type="checkbox"/> AIRIS Vento	Hitachi Medi...	MRI	AIRIS Vento	V1.0
<input type="checkbox"/> AIRIS Vento LT	Hitachi Medi...	MRI	AIRIS Vento LT	V1.0


Add Close

- 4 Click  button next to Package File Name.




Select file to upload at Select File. Input File Display Name and Description. And then, click Upload button. Don't touch Permanent File check box.

- 5 Click Add button on Software Package dialog.




- 6 Select the device that restores and click [SERVICE] - [FILE MANAGEMENT] - [Transfer files].



- 7 Select "System files" as [From] [Browse].



8 Tick desirable service back up file to transfer to the device.




9 Specify file saving directory.

Enter "/Download "in [Target directory path:].



10 Click [OK].

- 11 Click [Refresh device files] after a while and select “Download” folder on File Transfer screen, so that transferred Service Backup .zip file can be confirmed.




NOTE Take five minutes before Refresh device files.


(2) Restore a backup file

Service backup that has been transferred (downloaded) to the device is restored. Operate according to the following procedure.

- 1 Select the device for Service restore and click [Run command].



- 2 Select [Service Restore] and click [Run].



- Open [COMMAND RESULTS] tab and executed command “Status” is “Started”.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
ImportSB	Service Restore	2006/09/1...		Started	hirota

- Close [COMMANDS] screen and click [Refresh] upper-right corner of the screen.

NOTE Take five minutes before updating.

- Click [Run command] again and open [COMMAND RESULTS] tab.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
ImportSB	Service Restore	2006/09/1...	2006/09/1...	Completed	hirota
GetSB	Service Backup	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Timed out	hirota
CaptScrn	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota
		2006/09/1...	2006/09/1...	Completed	hirota

INPUT PARAMETERS		OUTPUT PARAMETERS	
Parameter	Value	Parameter	Value
		RtnImportSB	ok

- 6 “Completed” character is displayed in “Status” when process is completed. Select run commands.

PARAMETERS		COMMAND RESULTS			
Command	Description	Start Date	End Date	Status	By
ImportSB	Service Restore	2006/09/1...	2006/09/1...	Completed	hirota


- 7 “ok” is displayed in the value of [OUTPUT PARAMETER] when process is normally completed.

OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	ok

- 8 In order to make the contents of the service backup available, please turn off the MRI system and turn it on again.

7.4.2.5 Disk Cleanup


Alerts occur if device BE PC/ Free Space C: or FE PC/ Free Space C: is under threshold. Free Space needs to be increased.

<input type="checkbox"/> BE PC/HDD Size C:	 73402363904 Byte
<input type="checkbox"/> BE PC/HDD Free Space C:	 1066384128 Byte
<input type="checkbox"/> BE PC/HDD Life Time	 1982 Hour


NOTE Not only BE PC/ Free Space C: but also FE PC/ Free Space C:, FE PC/ Free Space D: are performed all at once.

Operate according to the following procedure.

- 1 Select the device to perform Disk Cleanup. Click [Run command].



2 Select [Disk Cleanup] and click [Run].



3 Open [COMMAND RESULTS] tab and executed command "Status" is "Started".

Command	Description	Start Date	End Date	Status	By
RemoteCle...	Disk Cleanup	2006/09/1...		Started	hirota

4 Close the command screen once, take 5 minutes and then click [Refresh] button at the upper right of the screen.

5 Click [Run command] again and open [COMMAND RESULTS] tab. "Completed" character is displayed in "Status" when process is completed. Select this item because abnormal or normal can not be assessed at this point.

Command	Description	Start Date	End Date	Status	By
RemoteCle...	Disk Cleanup	2006/09/1...	2006/09/1...	Completed	hirota

- 6 “ok” is displayed in the value of [OUTPUT PARAMETER].

OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	ok

- 7 Free space is increased when process is normally completed.


<input type="checkbox"/> BE PC/HDD Size C:	 73402363904 Byte
<input type="checkbox"/> BE PC/HDD Free Space C:	 60287528960 Byte
<input type="checkbox"/> BE PC/HDD Life Time	 1982 Hour

NOTE Not only BE PC/ Free Space C: but also FE PC/ Free Space C:, FE PC/ Free Space D: are performed all at once.


7.4.2.6 Shutdown

If some sorts of system problems occur and system should be shot down, system shot down is possible.

- 1 Select the device to perform Shutdown and click [Run command].




- 2 Choose [Shutdown] and enter password.



NOTE Inform Remote Service Administrators of "Unique Device Name" (MRI_ECHELON_Smart_Serial number) in advance and inquire passwords. Issued passwords are valid within 24 hours.

3 After entering password, click [Run] button.



4 By opening [COMMAND RESULTS] tab, "Status" of the executed command becomes "Started".

Command	Description	Start Date	End Date	Status	By
Shutdown	Shutdown	2006/10/1...		Started	hirota

5 Close the command screen once, take 1 or 2 minutes and then click [Refresh] button at the upper right of the screen.

6 The device is powered off after [Shutdown] was normally completed.

- 7 By choosing the executed command from [COMMAND RESULTS] tab, [OUTPUT PARAMETERS] value is “ok”.

OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	ok

NOTE [OUTPUT PARAMETERS] value is “Invalid Password” if wrong password was entered when a command was executed.

OUTPUT PARAMETERS	
Parameter	Value
Return result of ...	Invalid Password

7.4.2.7 Collect Manage Support Log


Manage Support Log can be acquired.

NOTE If has not been measured since started the device can not be acquired.


Acquire Manage Support Log according to the following procedure.

(1) Collect Manage Support Log execution

- 1 Select the device for Manage Support Log and click [Run command].



- 2 Select [Collect Manage Support Log] and click [Run] button.




- 3 The Manage Support Log file “ManageSupportTotal.xml”, “ManageSupportDaily_{{Year-Month-Day-Hour-Minute-Second}}.xml” is automatically uploaded on the server after a few minutes.


(2) How to download and display from server

If downloading Manage Support Log from server to user's PC, operate as follows.


- 1 Select the device where Manage Support Log was carried out and click [SERVICE] - [FILE MANAGEMENT] - [Manage system files].




- The screen of System Files is displayed, so select “ManageSupportTotal.xml” or “ManageSupportDaily_{"Year-Month-Day-Hour-Minute-Second"}.xml” and click [Download] button.




- The following window is displayed, so click a link.




- 4 The dialog of download displayed, so click [Save] button.



- 5 Specify saving place, the file name and click [Save] button.



- 6 [Download Complete] dialog is displayed, so click [Close] button.




7.4.2.8 Get Prescan Data

Prescan Data can be acquired.


Acquire Prescan Data according to the following procedure.

(1) Prescan Data execution

- 1 Select the device for Prescan Data and click [Run command].



- 2 Select [Get Prescan Data] and click [Run] button.




- 3 The Prescan Data file
“PrescanData_{"Year-Month-Day-Hour-Minute-Second"}.zip” is automatically uploaded on the server after a few minutes.


(2) How to download and display from server

If downloading Prescan Data from server to user's PC, operate as follows.


- 1 Select the device where Prescan Data was carried out and click [SERVICE] - [FILE MANAGEMENT] - [Manage system files].



- 2 The screen of System Files is displayed, so select "PrescanData_{\"Year-Month-Day-Hour-Minute-Second\"}.zip" and click [Download] button.




- 3 The following window is displayed, so click a link.




- 4 The dialog of download displayed, so click [Save] button.



- 5 Specify saving place, the file name and click [Save] button.



- 6 [Download Complete] dialog is displayed, so click [Close] button.




7.4.2.9 Password Reset


Change the password to the password of that year is possible if accidentally changed the date the following year.

Password reset according to the following procedure.

- 1 Select the device for change Password and click [Run command].




- 2 Choose [Password Reset] and enter password.



NOTE Inform Remote Service Administrators of "Unique Device Name"

- 3 After entering password, click [Run] button.



- 4 System restart after a few minutes.
- 5 Enter the password for that year at the login screen and log on.


7.4.2.10 Collect Hardware Status Log

To investigate device at the time of system down and failure, hardware status log can be acquired.


Acquire system log according to the following procedure.

(1) Collect Hardware Status Log execution

- 1 Select the device for Collect Hardware Status Log and click [Run command].



- 2 Select [Collect Hardware Status Log] and click [Run].



- 3 Open [COMMAND RESULTS] tab and executed command "Status" is "Started".

Command	Description	Start Date	End Date	Status	By
GetHardwareSt...	Collect Hardware...	2006/09/1...		Started	hirota

- 4 Close COMMANDS screen and click [Refresh] upper-right corner of the screen.

NOTE Take five minutes before updating.

- 5 Click [Run command] again and open [COMMAND RESULTS] tab.

The screenshot shows the 'COMMAND RESULTS' tab with a table of command execution results. Below the table are two sub-tables: 'INPUT PARAMETERS' and 'OUTPUT PARAMETERS'.

Command	Description	Start Date	End Date	Status	By
GetHardwareSt...	Collect Hardware...	2006/09/1...	2006/09/1...	Completed	hirota
GetLog	Collect Log	2006/09/1...	2006/09/1...	Timed out	hirota
CaptScrn	Screen Capture	2006/09/1...	2006/09/1...	Completed	hirota
		2006/09/1...	2006/09/1...	Completed	hirota

INPUT PARAMETERS	
Parameter	Value

OUTPUT PARAMETERS	
Parameter	Value
RtnGetHard...	ok

- 6 “Completed” character is displayed in “Status” when process is completed. Select run commands.

The screenshot shows the 'COMMAND RESULTS' tab with a table of command execution results. The 'Completed' status and the 'GetHardwareSt...' command are highlighted with red boxes.

Command	Description	Start Date	End Date	Status	By
GetHardwareSt...	Collect Hardware...	2006/09/1...	2006/09/1...	Completed	hirota

- 7 “ok” is displayed in the value of [OUTPUT PARAMETERS] when process is normally completed.

The screenshot shows the 'OUTPUT PARAMETERS' table with the 'ok' value highlighted with a red box.


OUTPUT PARAMETERS	
Parameter	Value
RtnGetHard...	ok

- 8 The log file named {Unit name}Log_{"Year-Month-Day_Hour-Minute-Second"}.zip is automatically uploaded into the server in a few minutes.


(2) Way of downloading from server and displaying

Download Hardware Status Log to pc you use from server at the time of log file reference.
Operate according to the following procedure.


- 1 Select the device for Collect Hardware Status Log and click [SERVICE] - [FILE MANAGEMENT] - [Manage system files].




- 2 The screen of System Files is displayed, so select {Unit name}Log_{"Year-Month-Day_Hour-Minute-Second"}.zip and click "Download".




- 3 The following window is displayed, so click a link.



- 4 The dialog of file download is displayed, so click “Save”.



- 5 Specify saving place, the file name and click [Save].



- 6 Execute the same operation to download other status log files.

7.5 Add readings

[Add readings] can be used to enter default at the registration of the device to the server or to manually enter the value in property.

It is possible to manually enter the value into monitoring database that compiles obtained data.

The image shows a software interface for adding readings. On the left is a sidebar menu with the following sections:

- SERVICE**
- ACTIONS**
 - Manage alerts
 - Manage schedules
 - Run command
 - Add readings** (highlighted with a red box)
- TOTAL ACCESS**
 - Connect to device
- FILE MANAGEMENT**
 - Refresh device files
 - Transfer files
 - Manage system files
 - View transfer log

An arrow points from the 'Add readings' option to the 'Add Readings' dialog box on the right. The dialog box has a title bar with a question mark icon. It contains a table with two columns: 'Property' and 'Value'. The table has five rows:

Property	Value
SL000000_Helium_Alarm	<input type="text"/>
SL000001_Helium_War...	<input type="text"/>
SL000002_EIS_Drive_...	<input type="text"/>
SL000003_ERDU_Batte...	<input type="text"/>
SL000004_ERDU_Batte...	<input type="text"/>


Below the table is a text area labeled 'Reason for entry'. At the bottom right of the dialog box are two buttons: 'Send' and 'Cancel'.

Figure7-40 Add readings

Monitoring object items are listed.




Figure7-41 Add readings

- (1) Property
Monitoring target property name is displayed.
- (2) Value
Enter the value to add.
- (3) Reason for entry
The reason for Add readings is entered.
- (4) Send
The added value is sent.
- (5) Cancel
Add readings is canceled and the screen is closed.
- (6) Help
Click  , so online help is displayed.

The operation procedure is described below.

- 1 Enter change value in desirable property input field.



SL000400_Helium_Level 100

- 2 Enter the reason for entry.



Reason for entry Change Value

- 3 Click [Send] button.



Send Cancel

- 4 Changes are reflected in Reading screen.

<input type="checkbox"/> MCU/Communication	OK
<input checked="" type="checkbox"/> MCU/Helium Level	100.0 %
<input type="checkbox"/> MCU/Magnet Internal Pressure	15.50 psiA
<input type="checkbox"/> MCU/Heater Current	0.20 W

7.6 File Management


7.6.1 Refresh device files

Use [SERVICE] - [FILE MANAGEMENT] [Refresh device files] if confirming the information of files which exist in the system as follows.


- To transfer files from the device to the server, confirmation of files which exist in the device.
- After they were transferred from the device to the server, confirmation of that files were transferred.

Operate according to the following procedure.

- 1 Click [SERVICE] - [FILE MANAGEMENT] [Refresh device files].



- 2 "Device File Update" of the following screen is displayed.



Device File Update

DEVICE FILE FILTER SCHEDULE

Device MRI_ECHELON_Smart_40917


From directory /

Filter *

Subdirectories Include

Schedule Cancel

- 3 Click "Schedule".



NOTE Use default. Do not change tabs and change setting.

- 4 Minutes later, files of the system can be referred to from [SERVICE] - [FILE MANAGEMENT] [Transfer files] .

7.6.2 Transfer files

[Transfer files] can be used to transfer from the server to the device or from the device to the server.

Select the device to transfer files and click [SERVICE]-[FILE MANAGEMENT]-[Transfer files].




Figure7-42 Transfer files

7.6.2.1 Transfer files from the device to the server




Figure7-43 Screen for transfer from the device to the server

Table7-11 Setup Items of File Transfer (MRI device to server)


1. Screen switching	TRANSFER screen and SCHEDULE screen can be switched.
2. From	Transfer source is checked.
3. Directory	Transferer directory is displayed.
4. File	A file in the selected directory is displayed.
5. Last Modified	The last modified date is displayed.
6. Size	File size is displayed.
7. Send	A file is transferred.
8. Cancel	Transfer is canceled and this screen is closed.
9. Delete files after upload	Tick here if deleting a file after uploading.

The following is the procedure to transfer from the device to the server.


- 1 Tick “Device” Browse box on File Transfer screen.




- 2 Click a directory that includes desirable file to transfer to the server.



3 Tick desirable file to transfer to the server.



4 Click [Send].



NOTE If you delete a file in the device after transferring a file to the server, tick Delete files after upload before transferring.

7.6.2.2 Transfer files from the server to the device




Figure7-44 Screen for transfer from server to the device

Table7-12 File transfer setup items (the server - the device)


1.	Screen switching	TRANSFER screen and SCHEDULE screen can be switched.
2.	From	Transferer is ticked.
3.	File	The forwarding file is ticked.
4.	Description	File description is displayed.
5.	Size	File size is displayed.
6.	Last Modified	The last modified date is displayed.
7.	To	It is automatically changed after checking a transferer.
8.	Directory	Transferer directory is displayed.
9.	File	A file in the selected directory is displayed.
10.	Last Modified	The last modified date of file is displayed.
11.	Size	File size is displayed.
12.	Target directory path	Specify the directory that a file is saved in. (Specification is possible by directly clicking an icon in a folder)
13.	Overwrite files	Choose whether to overwrite if same name file already exists.
14.	Create directory	New directory is created.
15.	Send	File is transferred.
16.	Cancel	The transfer is canceled and this screen is closed.

The following is the procedure to transfer from the sever to the device is described below.

- 1 Tick “System files” Browse box of [From].




- 2 Tick desirable file to send to the device.



NOTE Upload the file in advance in reference to section 7.6.3.1 if transferring files in Local PC.

3 Specify a directory that a file is saved in.

Enter “/Download “in [Target directory path] box.



4 Click [Send].

NOTE If the file that has same name already exists, tick Overwrite files before transferring.

7.6.3 Manage system files

The list of files stored in the server is displayed.

It is also possible to upload files in Local PC and download files stored in the server.

System Files in the following list are displayed.

Table7-13 Display of System Files

Condition	Display of Files
A file that is uploaded from user’s Local PC.	Yes
A file that is transferred from the device registered with the group user belongs to.	Yes
A file that is uploaded from Local PC by the other user.	No
A file that is transferred from the device registered with the group user does not belong to.	No

Moreover, the system file has the file type that exists in the following table.

Table7-14 System File kinds

File Name	Contents
Screen_{Device Name_YMDHMS}.png	This file is uploaded after image capture command is executed. See section 7.4.2.1
diag_dump_info_{System name}_{Serial number}_{"MonthDay"}_{"HourMinuteSecond"}.zip	This file is uploaded after log collection command is executed. See section 7.4.2.2
BELogFile.zip	
ScanLogFile.Zip	
ReconLogFile.zip	
ServiceBackup.zip	This file is uploaded after service backup command is executed. See section 7.4.2.3
DICOM{YMDHMS}.zip	This file is generated when the library is acquired from MRI device. See section 7.7.6
{User + “Serial No” + “Data”}.zip	{User + “Serial No” + “Date”}.zip is generated when a library is acquired from MRI device. See section 7.7.7
{YMD}.txt	{YMD}.txt is the defect report sent from MRI device by Sentinel Messenger. See Instruction Manual.
LogReader.msi	LogReader.msi is the install file of Log Reader.
{Manual Name}.pdf	{Manual Name}.pdf are various manual.
PrescanData_{"YMDHMS"}.zip	This file is uploaded after Get Prescan Data command is executed. See section 7.4.2.8

File Name	Contents
{Unit name}Log_{"YMDHMS"}.zip	This file is uploaded after Collect Hardware Status Log command is executed. See section 7.4.2.10
ManageSupportTotal.xml	This file contains operation report for one month. This file is generated once a month.
ManageSupportDaily_{"YMDHMS"}.xml	This file contains last operation report. It is generated every time power on the MRI system.

Select the device to refer files and click [SERVICE]-[FILE MANAGEMENT]-[Manage system files].




Figure7-45 Whole screen of Manage system files

The system file management screen is displayed.







Figure7-46 Manage system files (Device)

- (1) File
File names are indicated. Tick the box on the left side if downloading or deleting files.
- (2) Owner
File owner is displayed.
- (3) Size
File size is displayed.
- (4) Last Modified
The file date is displayed.
- (5) MD5 Checksum in Hex
File MD5 checksum is indicated.
- (6) Change display
File upload source which "System Files" screen indicates can be changed by "My Files" and "Device Uploaded Files".

- (7) Add
File can be uploaded from local PC if you click .
- (8) Delete
A ticked file is deleted if you click .
- (9) Download
By clicking, selected files can be downloaded.
- (10) Cancel
This operation is canceled and the screen closes.
- (11) Help
Online help is displayed if you click .

By changing displays, the list of files uploaded by users is displayed.







Figure7-47 Manage system files (User)


- (1) File
File names are indicated. Tick the box on the left side if downloading or deleting files.
- (2) Description
The description which was input at the time of upload is indicated.
- (3) For
The names of the device and the group which were selected at the time of uploading files are indicated. "Friendly name" which was registered in "Personal Configuration" is indicated at the time of upload without selection. (Refer to Chapter 7.13.)
- (4) Size
File size is displayed.
- (5) Last Modified
The file date is displayed.

- (6) Change display
File upload source which “System Files” screen indicates can be changed by “My Files” and “Device Uploaded Files”.
- (7) Add
File can be uploaded from local PC if you click .
- (8) Delete
A ticked file is deleted if you click .
- (9) Download
By clicking, selected files can be downloaded.
- (10) Cancel
This operation is canceled and the screen closes.
- (11) Help
Online help is displayed if you click .


7.6.3.1 Upload from Local PC to the server

The following is the procedure to upload from Local PC to the server.


- 1 Display “My Files” and click .




- 2 The reference screen of files to upload is displayed.



3 Click **Browse...**.




4 Select a file.




File name corresponds to only alphabetic one byte characters. Japanese file name can not be specified.

5 Input "Description" and click **Upload**.




6 The file is uploaded.




7.6.3.2 Download files from the server to Local PC.

The following is the procedure to download files from the server to Local PC.

- 1 Select the desirable file to download and click “Download”. (Multiple selection is accepted.)




- 2 The following window is displayed, so click a link.




- 3 The dialog of file download is displayed, so click “Save”.



- 4 Specify saving place, the file name and click [Save].





- 5 [Download complete] dialog is displayed, so click [Close] button.




7.6.3.3 Server file deletion

Operate according to the following procedure if deleting the file which was uploaded to the server.

- 1 Tick the box of the file to delete and click .



- 2 The dialog of “Confirmation” is displayed, so click “Yes”.



- 3 The file is deleted.

7.6.4 View transfer log

It is possible to refer to the transfer log from the server to the device and from the device to the server.

Selected device View transfer log is displayed.




Figure7-48 Whole screen of View transfer log

(1) Action	(2) Status	(3) By	(4) Notes	(5) Date
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:28:45
FileTransfer	FileUploadCompleted	MRI_ECHELON_Smar...	Member with member...	2016/11/29 16:28:45
FileTransfer	FileUploadStarted	MRI_ECHELON_Smar...	documentName=/U...	2016/11/29 16:28:45
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:54
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:54
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:54
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:53
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:52
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:52
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:51
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:51
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:50
FileTransfer	FileUploadCompleted	MRI_ECHELON_Smar...	Member with member...	2016/11/29 16:26:50
FileTransfer	FileUploadStarted	MRI_ECHELON_Smar...	documentName=/U...	2016/11/29 16:26:50
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:50
FileTransfer	FileUploadCompleted	MRI_ECHELON_Smar...	Member with member...	2016/11/29 16:26:49
FileTransfer	FileUploadStarted	MRI_ECHELON_Smar...	documentName=/U...	2016/11/29 16:26:49
FileTransfer	FileTransferCompleted	MRI_ECHELON_Smar...	Received upload file...	2016/11/29 16:26:49
FileTransfer	FileUploadCompleted	MRI_ECHELON_Smar...	Member with member...	2016/11/29 16:26:49
FileTransfer	FileUploadStarted	MRI_ECHELON_Smar...	documentName=/U...	2016/11/29 16:26:49

Figure7-49 View transfer log

- (1) Action
The action content is displayed.
- (2) Status
The action status name is displayed.
- (3) By
The transfer source is displayed.
- (4) Notes
Transfer Log is displayed.
- (5) Date
File transfer date is displayed.

7.7 Connect to device

- By remote desktop function, device can be controlled by remote and operation procedure can be explained.
- Because operation screen can be confirmed at the time of error occurrence, trouble situation can be grasped immediately.
- It is possible to confirm the operation screen as well as the taking picture parameter when a defective image is created.

Perform it according to the following flow if performing remote desktop in devices.





Figure7-50 A flowchart to perform remote desktop connection

7.7.1 Before Start


- 1 Click [SERVICE] - [TOTALACCESS] - [Connect to device].
The ActiveX install screen is displayed for Windows 7 before the remote desktop function is executed. (Initial connection only)
- 2 Select [Install ActiveX Control].



- 3 Click [Install] button.



4 Device can be connected if installation is completed.




7.7.2 How to start


Remote desktop is used on the device.

The following is the procedure.


- 1 Click [SERVICE] - [TOTAL ACCESS]-[Connect to device].




- 2 Pop-up window is displayed.




3 Wait until “Ready” character is indicated.




4 Click [Connect] button.



5 The following screen is displayed if proxy server is set up. Click [OK] button.




NOTE If “Automatically detect settings” is ticked in Internet Explorer LAN settings , the connection is unavailable. Set up the address and port.



- 6 “Security Alert” pop-up is displayed if Firewall in Windows 7 is available. Click “Unblock”. (Initial connection only)




- 7 Click [OK] if the next pop-up appears.



- 8 The following confirmation screen is displayed on the screen of the device.

Access is available only if a user clicks [Accept]. (Reject is automatically applied if a user does not click [Accept] in five minutes.)




A dialog that indicates the users connecting with devices is displayed if [Accept] button is clicked.




- 9 Login/Password screen is displayed if the connection is authorized by the device.

Enter a password. (Contact the system manager regarding a password.)




NOTE Close the screen and restart from the beginning if wrong password is inputted.

- 10 The viewer is launched when the connection is provided.
It takes few minutes to display the desktop.





11 Device desktop is displayed.



NOTE Wallpaper becomes black at the time of remote desktop connection.

12 Click or icon to terminate remote desktop connection.

NOTE If the following error appears at VNC connection, reconnect after a few minutes.



7.7.3 VNC viewer.....

The following VNC viewer is launched when the remote desktop is executed.




Figure7-51 whole screen of VNC Viewer

Main function tool bar is displayed at the top of the Viewer.






Figure7-52 Icons

Operate the following procedure if the screen size and color are changed.

- 1 Click an icon.




- 2 The pulldown menu is displayed.





- 3 Screen design can be changed by the pulldown menu.
The default screen size is “Auto Scale screen” and the default color is “Full Colors”.

- 4 An image is displayed at the same magnification if [Normal screen] is selected. An image is displayed on 256 color display if [256 Colors] is selected.



- 5 Screen design is changed.



Regarding VNC viewer, if line speed becomes slow, it sometimes causes a delay in the refresh of screen information. Refresh a screen manually if a screen is distorted. Click Refresh button  of VNC tool bar if refreshing a screen manually.




Figure7-53 Refresh Screen

7.7.4 Personal information protection mode.....

All items related to the following personal information are hidden so that personal information leak is prevented at the time of remote desktop connection.

Table7-15 Personal Information Protection Items

No	Personal information	How to hide
1	PatientID	Hidden with **** The latest Study ID is indicated as alternatives on Patient Directory screen.
2	AccessionNumber	Hidden with ****
3	AdditionalPatientHistory	Hidden with ****
4	BranchOfService	Hidden with ****
5	ContrastAllergies	Hidden with ****
6	CountryOfResidence	Hidden with ****
7	EthnicGroup	Hidden with ****
8	IssuerOfPatientID	Hidden with ****
9	MedicalAlerts	Hidden with ****
10	MedicalRecordLocator	Hidden with ****
11	MilitaryRank	Hidden with ****
12	OtherPatientIDs	Hidden with ****
13	PatientAddress	Hidden with ****
14	PatientAge	Hidden with ****
15	PatientBirthDateTime	0001/1/1 00: 00: 00.000
16	PatientComments	Hidden with ****
17	PatientInstitutionResidence	Hidden with ****
18	PatientMotherBirthName	Hidden with ****
19	PatientName	Hidden with ****
20	PatientReligiousPreference	Hidden with ****
21	PatientSex	Other
22	PatientSize	-
23	PatientTelephoneNumbers	Hidden with ****
24	PatientWeight	-
25	PerformingPhysiciansNameList	Hidden with ****
26	PregnancyStatus	Unknown

No	Personal information	How to hide
27	Radiologist	Hidden with ****
28	RegionOfResidence	Hidden with ****
29	Smoking Status	Unknown
30	Technologist	Hidden with ****

The following operation accompanying inputting and referring to personal information is restricted.

- Correct Patient
- Patient Registration
- Patient Information Property (It is possible to display only.)

But it is possible to use “Rapid Registration”, register the patient and scan.

Hiding process target screens are listed below.

Table7-16 Screens that are subject to hiding process

No	Screen
1	Patient Directory
2	Patient Information Property
3	Patient Registration
4	Exam Window
5	Study Selector
6	Scan UI
7	Series Selector
8	Task View
9	Accounts Setting
10	TM Raw Data Selector
11	TM Prescan Data Log Selector
12	TM Series Selector
13	Patient Registration Setting

7.7.5 Service Tool & Diagnostic.....

It is possible to use “Service Tool & Diagnostic” by remote desktop connection. Operate the following procedure.

- 1 Click [Maintenance] and "Service Tool & Diagnostic" is displayed.




Figure7-54 Maintenance menu

- 2 Select and start tools from "Service Tool & Diagnostic" menu.

NOTE It is not necessary to log in with [User Recognition] because license is automatically registered during remote desk top connection.

7.7.6 How to acquire and display image.....

If there is something wrong with image, image can be acquired from device and it can be displayed.

(1) Image upload

The following is image upload procedure.

- 1 Click [Patient Directory].



Figure7-55 Patient directory menu

- 2 The following screen is displayed.




Figure7-56 Patient directory screen

- 3 Select the patient to upload from Patient List and choose [Remote Server] from [Send] of the menu displayed by right-click.




Figure7-57 Uploading patient image data

- 4 Click [OK] to transfer image data to support center server.




Figure7-58 File Upload


(2) Image download

The image uploaded to server is uploaded into local PC. Operate as follows.


- 1 Choose the device where image has uploaded, and select [SERVICE] - [FILE MANAGEMENT]-[Manage system files].




- 2 The screen of System Files is displayed, so select “DICOM{YearMonthDayHourMinuteSecond}.zip” and click “Download”.




- 3 The following window is displayed, so click a link.




- 4 The dialog of file download is displayed, so click "Save".



- 5 Specify saving place, the file name and click [Save].



- 6 [Download Complete] dialog is displayed, so click [Close].



(4) Image display

It is possible to extract downloaded image and display it with attached viewer software. Operate according to the following procedure.

- 1 Open explorer, confirm the folder where the file downloaded in (3) has been saved and file is being saved.




Figure7-59 My documents folder

- 2 Right-click screen.




Figure7-60 My document folder select

3 Select and click [New] - [Folder].




Figure7-61 Folder creation

4 "New folder" is created, so change its name to any name.




Figure7-62 Folder name input

5 Use ZIP file extraction tool and copy file to the new folder.




Figure7-63 File extraction

6 Double-click [VIEWER].




Figure7-64 VIEWER folder

7 Double-click [DCMVIEW.exe].




Figure7-65 Execute DCMVIEW.exe

8 License confirmation screen is displayed. Click [OK] button if agreeing with description.



Figure7-66 Authentication screen

9 Image is displayed as follows.




Figure7-67 Start screen

See Help for detailed usage.

7.7.7 How to get library.....

It is possible to get library from device.

(1) Library upload (Device to Server)

Library is uploaded once from device to server for getting library. Operate the following procedure.

- 1 Click [Protocol/Task Library].



Figure7-68 Menu protocol/Task library

- 2 The following screen is displayed.




Figure7-69 Library tree screen

- 4 Right-click and select [Send]-[RemoteServer].




Figure7-71 Library export

- 5 Compression file [User + “Serial No” + “Date”.zip] is automatically uploaded to support center server.






Figure7-72 Upload completion


- 6 Select the device where Library upload was carried out and click [SERVICE] - [FILE MANAGEMENT]-[Manage system files].



- 7 The screen of System Files is displayed, so select [User + “Serial No” + “Date”.zip] and click “Download”.



- 8 The following window is displayed, so click a link.




- 9 The dialog of file download is displayed, so click "Save".




10 Specify saving place and click [Save].

NOTE Do not change the name of the file. The file can not be imported if the name was changed.




11 [Download complete] dialog is displayed, so click [Close] button.




By using a viewer, a downloaded protocol/task library can be referred to. According to the following procedure, refer to it.


1 Extract a downloaded file and double click “AutoRun.exe”.



2 A viewer boots up.




NOTE If Microsoft .NET Framework Version 1.1 is not installed, the following error dialog box appears. Install Microsoft .NET Framework Version 1.1.




(2) Import to a system of Support Center

Regarding a library which was acquired from a system of customer site, you can import it to a system of Support Center and investigate it. According to the following procedure, operate.


- 1 Into a CD-R / DVD, save a library which was downloaded from a server.




- 2 Into the system, insert the CD-R or DVD where the library has been saved.
- 3 By clicking [Maintenance], "Service Tool & Diagnostic" appears.




- 4 "System Control" tab becomes active. So choose "Protocol Import" from a pull-down menu.




- 5 "Protocol File Import" screen appears.




NOTE Unless the CD-R or DVD, which the library has been saved in, is inserted, nothing is indicated.




- 6 Choose a library to be imported, and click "Import".




- 7 The library is imported.



- 8 By clicking [File] and [Exit], close "Protocol File Import" screen.



- 9 On the [Protocol/Task Library], the library that was mounted to "Remote Server" is imported. Right-click [Send] and select "LocalDB".



NOTE Transfer to "LocalDB" before the system is powered off. If the system is powered off, the library disappears from "Remote Server".


7.8 Update software

Service pack is applicable to device.

System administrator uploads service pack into server.

7.8.1 Displaying [Software update] Dialog.....


- 1 Select Service pack applicable device.



- 2 Click [SERVICE] - [Sentinel Update] - [Update software].




3 Open [Software Updates].




7.8.2 Service pack applicable schedule

Select applicable service pack at first and register application schedule.


- 1 Select applicable service pack from [Software Updates] - [SOFTWARE PACKAGES].




- 2 Lower side of the dialog indicates service pack information. See the information indicated in [ASSIGNED DEVICE TYPES] tab and [PROFILE] tab, and make sure indicated service pack is applicable.




- 3 Click [SCHEDULE] tab.



- 4 Create new schedule.
Click  button.



- 5 [New Schedule] dialog opens.
Enter easy name into [Name] and click [Add] button.




NOTE If name is registered once, it can not be used twice. Use different name by adding an identifier such as a number to the file name and so on, even if deleting schedule and hoping to use same name again.

Reference: Item description

Name	Only this is mandatory item. See above NOTE . Up to 50 characters (25 two-byte characters).
Description	Schedule description. Up to 200 characters (100 two-byte characters).
Type	Device task and Enterprise task are selectable. Device task is device initiative, so center load is low, but it has something hard to manage. Enterprise task is center initiative, so it is easy to manage, but center sever load is high and delay often occurs.
Frequency	It is possible to repeat schedule, but make service pack application frequency once in general.
Start	It is possible to set up application start date. NOTE: Service pack is not always applied on application start date due to device status.
Remind	If service pack is not applied up to the date set here, alarm occurs and if mail notification is set up, mail is received.
Install by	As the case of Remind, alarm occurs if service pack is not applied up to the date set here and if mail notification is set up, mail is received. Separate to use them, like setting up “Remind” before application time limit and setting up “Install by” after application time limit.

6 Schedule is registered.



7.8.3 Application result confirmation

The service pack is forwarded to the device, and installed automatically when becoming a date when the schedule registration was done.

Service pack application result can be confirmed with [Installed Software].

MRI_ECHELON_Smart_40925 Assigned to ▾

● MRI_ECHELON_Smart_40925 Matched configuration

ECHELON Smart V6.0D

Actual configuration information Nov 29, 2016 4:25:44 PM

Name	Value
Type 5	Type 5
V6.0D	V6.0D
DWI Package	DWI Package
BASG Imaging	BASG Imaging
Water Excitation	Water Excitation
RAPID	RAPID
Service Pack (ECHELON_Smart)	Sentinel Certification 2013 SP

When the installation succeeds, the service pack name is displayed.

Figure7-73 Application result confirmation

That's not to say that "Installed Software" information is updated as soon as [Software Updates] operation described in chapter 7.8. It is regularly updated, but perform [Get actual configurations], (see chapter 7.9), if you hope to acquire information immediately.

7.9 Get actual configurations

Device "Installed Software" information is regularly updated. But by using [SERVICE] - [Sentinel Update] - [Get actual configurations], device "Installed Software" information can be updated.






Figure7-74 Device summary

The following is "Get actual configurations" usage procedure.

- 1 Select the device acquiring "Installed Software" information and click [SERVICE] - [Sentinel Update]-[Get actual configurations].




- 2 “Actual Configurations” screen is displayed, so click .



Task	Type	Start Date	Status
Collect revision for ...	Enterprise task	2006/08/01 15:54:01	Active
Collect revision for ...	Enterprise task	2006/08/01 15:54:39	Stopped

Close

- 3 “New Schedule” screen is displayed.



Description

Name *

Description

Type

Recurrence

Frequency


Start

Start / / :


Add Cancel

4 Fill in [Name], choose [Type] and click [Add].

NOTE Description omission is acceptable.



5 Make sure task has been added and then click [Close].



- 6 Take a minute to update and the latest "Installed Software" information is indicated.

The screenshot displays the configuration page for an MRI system. At the top, a header bar shows the system name 'MRI_ECHELON_Smart_40925' and a status 'Assigned to'. A callout box points to the text 'Date is updated.' which is positioned above a red-bordered box containing the date and time 'Nov 29, 2016 4:25:44 PM'. Below this, a table lists the 'Actual configuration information' for the system.

Name	Value
Type 5	Type 5
V6.0D	V6.0D
DWI Package	DWI Package
BASG Imaging	BASG Imaging
Water Excitation	Water Excitation
RAPID	RAPID
Service Pack (ECHELON_Smart)	Sentinel Certification 2013 SP

7.10 View chart

By using [SERVICE] - [REPORTS] - [View chart], acquired data can be listed and graphed on screen.

The following is the procedure up to View chart usage.

1 Select device.

The screenshot displays the Sentinel Customer Support web application in a Microsoft Internet Explorer browser. The page title is "Sentinel Customer Support" and the user is signed in as "Hitachi". The breadcrumb navigation shows "All > ECHELON Smart Devices > MRI_ECHELON_Smart_40925".

The main content area shows the configuration for the device "MRI_ECHELON_Smart_40925". It includes a "Matched configuration" section with a table of parameters:

Parameter	Value
Name	Type 5
Type 5	V6.0D
V6.0D	DWI Package
DWI Package	BASG Imaging
BASG Imaging	Water Excitation
Water Excitation	RAPID
RAPID	CE-MRA Package
CE-MRA Package	

Below the configuration is a "Readings" section with a table of 299 properties. The table includes columns for the property name, a bar chart, a numerical value, and a timestamp. The "Readings" section is currently empty, with a "Select All" checkbox and a "299 Properties" label.

The right-hand sidebar contains several sections: "SERVICE" (with actions like "Manage alerts", "Manage schedules", "Run command", "Add readings"), "TOTAL ACCESS" (with "Connect to device"), "FILE MANAGEMENT" (with "Refresh device files", "Transfer files", "Manage system files", "View transfer log"), "SENTINEL UPDATE" (with "Update software", "Get actual configurations"), and "REPORTS" (with "View chart", "View report").

The bottom of the page shows "Controls", "Activities" (with "Last activity: FileRepositoryManager uploadFileNotify (Nov 29, 2016...)"), and "Alerts" (with "0 Alerts").

- 2 Click [SERVICE] - [REPORTS] - [View chart].



- 3 Select [SERVICE] - [REPORTS]-[View chart] and the following screen is displayed.








Figure7-75 View chart

- (1) Properties
Select the property to indicate. (Selection is mandatory. Multiple selections are possible. Selecting up to 6 items is possible.)
- (2) Devices for comparison
Select comparing devices. (Selection is mandatory. Multiple selections are possible. Selecting up to 6 items is possible.)
- (3) Filter
Filtered and displayed value can be set.
- (4) Maximum number of readings
Displayed maximum number is entered. (Default is 1000)
- (5) Start
Start date is specified.
- (6) End
End date is specified.
- (7) Submit
By clicking , chart or table is displayed.
- (8) Display switch
By clicking , chart or table display is switched.
- (9) Close
By clicking , this screen closes.
- (10) Help
By clicking , on-line help is displayed.
- (11) Display screen
Chart or table is displayed here.

NOTE [Properties]×[Devices for comparison]≤6 is number of selectable items in [Properties] and [Devices for comparison].

The following is View chart procedure.

- 1 Tick items to display in [Properties] and [Devices for comparison].

NOTE Select same unit properties if choosing two or more and displaying chart.

<input type="checkbox"/>	BEPC/HDD Size D:
<input checked="" type="checkbox"/>	BEPC/HDD Free Space C:
<input type="checkbox"/>	BEPC/HDD Free Space D:
<input type="checkbox"/>	BEPC/HDD Life Time
<input type="checkbox"/>	BEPC/HDD Temperature

Devices for comparison

<input type="checkbox"/>	MRI_ECHELON_Smart_40917
<input checked="" type="checkbox"/>	MRI_ECHELON_Smart_40918
<input checked="" type="checkbox"/>	MRI_ECHELON_Smart_40919

- 2 Enter and select [Filter], [Maximum number of readings], [Start] and [End] if necessary.

Filter

MCU/Helium (Fault)	▼	=	▼	
--------------------	---	---	---	--

Maximum number of readings

1000

Start


8	/	24	/	2006			:	
---	---	----	---	------	--	--	---	--

End

8	/	31	/	2006			:	
---	---	----	---	------	--	--	---	--






Submit


3 Click and chart is displayed.




- (1) Properties
Graph display selection is possible.
- (2) Time period
By moving the cursor, data range specification is possible.
- (3) Graph display
X axis period zoom in is possible within selected range. Besides, graph display becomes normal one by double-clicking graph screen.
🕒 is the latest Alert of the device chosen in “Devices for comparison”.

4 Screen can be zoomed in as follows.



- 5 Click  and display is switched to table display.




- (1) Date
Acquisition time is indicated.
- (2) Value
Value is indicated.
Indicating items increases at the time of multiple selections.

7.11 Report


By acquired data, automatically generated report can be browses.

7.11.1 View report

- 1 If selecting any group, [SERVICE] – [REPORTS] - [View report] become active.



2 Click [View report].



3 Report menu is displayed.







Figure7-76 Report menu screen

- Click “Channels” of the menu on the upper left of the screen and display channel list.



- Select any channel from channel list.



NOTE “Personal Channel” is not available.

- 6 Channel inside displays the report list automatically generated and available for browse.





Figure7-77 Channels


- 7 Click name of report to browse and the report is displayed.

 [Device Status Information Report](#)

Device Status Information Report

 [Device Status Information Report](#)

Device Status Information Report



https://www.sentinel-service.com - Questra Smart Service - [Device Status Information Report] - Microsoft

TOC First Prev Next Last Goto Page 1 of 1 100% Search Download Print


HITACHI Device Status Informatic

Registered Devices

<i>Device Name</i>	<i>Model</i>	<i>Version</i>
MRI_ECHELON_Smart_40919	ECHELON Smart	V4.0A
MRI_ECHELON_Smart_40925	ECHELON Smart	V4.0A
MRI_ECHELON_Smart_40929	ECHELON Smart	V4.0A

Figure7-78 Report browse screen

Report operation is possible from the menu at the top of report browse.



- (1) TOC
Report content opens.
- (2) First
Page can turn back to first page.
- (3) Prev
Page can turn back to previous page.
- (4) Next
Next page opens.
- (5) Last
Last page opens.
- (6) Goto
Page entered in page opens.
- (7) Page
“Current page/Total page” is displayed.
- (8) Screen size
Screen size can be changed.
- (9) Search
Search bar is displayed.
- (10) Download
Report can be saved.
- (11) Print
Report is printed.
- (12) Help
On-line help is displayed.
NOTE Not available now.
- (13) Close
Report closes.

7.11.2 Using content

- 1 According to above procedure, display the report to browse.
- 2 Click “TOC” from the menu at the top of the report browsing screen and display content.








Figure7-79 Content display

- 3 Click the caption indicated on the left and its caption page can be displayed.

7.11.3 Using search function

- 1 According to above procedure, display the report to browse.
- 2 Click "Search" from the menu at the top of screen browsing report.



- 3 Screen left side displays search menu.





Figure7-80 Search menu display

- 4 Click the items to search from report body. In the following case, to search with device name, “MRI_ECHELON_Smart_40919” indicated at the top of Device Name items is clicked. By clicking an item, it is added to search menu on the left.




- 5 If entering any phrase and clicking [Search Now], report search starts according to the phrase.




- 6 By clicking indicated search result, jumping to the page including its phrase is possible.

- Click [New Search] if searching with other phrase again. Screen returns to the entry screen by clicking. New item can be added here, too. (Its way is the same as above item addition procedure.)




- Only administrator user can use [Save Search] button next to [New Search].
- By marking "*" to search phrase, all including phrase before * is indicated as search result. In the following case, to search device name, all devices whose name are same and whose serial numbers are different are searched.



7.11.4 Report save

Displayed report can be saved.

- 1 According to above procedure, display any report.
- 2 Click “Download” from the menu at the top of report browsing screen.



- 3 Download screen is displayed.




Figure7-81 Download screen

4 Saving format is selected.

5 type formats can be selected.

Export Report To:

PDF

PDF Quality: Split Large Pages:

Page Width: Page Height:

Excel Data

Excel Display

RTF

Fully Editable RTF

Tips:

1. Upgrade IE to version 4.0.5 (or later) or save the PDF report locally before viewing.
2. Excel Data format is good for data manipulation. It was designed for tabular and listing reports.
3. Fully Editable RTF format is good for multi-control editing, but creates significantly larger files than RTF format.
4. PDF Quality level 100 gives the lowest image quality but the smallest PDF file size, and 300 gives the highest image quality but the largest PDF file size.

(1) PDF

If selecting PDF, it is additionally possible to set up image quality, the process for the case where report does not fit into one page, page width and height. In addition, download is possible even without setting them up.

(2) Excel Data

Output file is Excel format, so its ruled line and so on sometimes looks a little bit different from real report.

(3) Excel Display

Output file is Excel format, but its ruled line and so on is the same as real report.

(4) RTF, Fully Editable RTF

If selecting RTF, rich text format file is output. It can be opened with Microsoft Word.

NOTE If a report is saved in “Excel data” and “Excel Display”, number of pages to download need to be 50 pages or less. If more than 50 pages are downloaded, “Warning: Excel converter page limit of 50 exceeded.” is indicated after 50 pages.

5 Set up downloading page range.

As the following figure description, any page can be specified, too.

Page Range:

All

Current page

Pages:

Enter page numbers and continuous page ranges separated by commas. For example: 1,3,5-12.


- 6 Click [View Report] or [Save Report]. If selecting [View Report], report immediately opens, but it is not saved. If selecting [Save Report], report does not open and is saved.



7.11.5 Report print.....

Displayed report can be printed.

- 1 According to above procedure, display any report.
- 2 Click "Print" from the menu at the top of report browsing screen.



- 3 Report is automatically converted to PDF format and print screen is displayed.




Figure7-82 PDF format report and print screen

- 4 Click [Print].

7.12 Personal Configuration

[Edit your profile] enables log-in user information display and change. Click [CONFIGURATION] - [PERSONAL CONFIGURATION] - [Edit your profile].




Figure7-83 Edit your profile





Figure7-84 User Profile

- (1) Display switch
[PROFILE] and [PROPERTIES] display is switched.
- (2) Save
Value changed setting is saved.
- (3) Cancel
This operation stops and this screen closes.
- (4) Help
By clicking , on-line help is displayed.
- (5) Profile
Profile is displayed.

User

(1) First name

(2) Last name

(3) Friendly Name

(4) Company

(5) Description

(6) User Name (7)

(8) Role

(9) Email address

(10) Phone number

Address

(11) Location

(12) Description

(13) Address1

(14) Address2

(15) City

(16) State or Province (17) Postal Code

(18) Country

Regional settings

(19) Time zone

User image

(20)

Figure7-85 Profile

- (1) First name
First name is indicated.
- (2) Last name
Last name is indicated.
- (3) Friendly name
Friendly name is indicated.
- (4) Company
Company is indicated.
- (5) Description
Department is indicated.
- (6) User Name
User Name is indicated. (NOTE: Administrator registers user.)


- (7) New password
Password can be changed.
- (8) Role
Role is indicated. (Administrator authority is required for Role setting change.)
- (9) Email address
Registered Email address is indicated.
- (10) Phone number
Registered phone number is indicated.
- (11) Location
Location is indicated.
- (12) Description
Registered address description is indicated.
- (13) Address1
Address1 of registered address is indicated.
- (14) Address2
Address2 of registered address is indicated.
- (15) City
Registered address city is indicated.
- (16) State or Province
Registered address state or province is indicated.
- (17) Postal Code
Registered address postal code is indicated.
- (18) Country
Registered address country is indicated.
NOTE Administrator user authority is required for setting change.
- (19) Time zone
Time zone is indicated.
- (20) User image
An image of one's face photograph and so on can be uploaded. (An image needs to be 120 pixels in width and 120 pixels in height, and its file format needs to be JPG.)

If clicking **Upload**, image of one's image can be uploaded. Follow the next procedure and upload it.


- 1 Click [CONFIGURATION] - [PERSONAL CONFIGURATION] - [Edit your profile].

- 2 Click **Upload**.

3 Click **Browse...**.




4 Choose a file.




An image needs to be 120 pixels in width and 120 pixels in height, and its file format needs to be JPG.

- 5 Click **Upload**.



- 6 The image is uploaded.



Click to change password.




Figure7-86 New password

- (1) Old password
Enter current password.
- (2) New password
Enter new password. (Password must be more than 5 characters and less than or equal to 64 characters).
- (3) Confirmed password
Enter input new password again.
- (4) OK
Password can be changed.
- (5) Cancel
This operation stops and this screen closes.

NOTE Errors occur if entry is wrong.

Property screen enables Extended Property setting.








Figure7-87 PROPERTIES

- (1) Add
By clicking , Extended Property can be added.
- (2) Edit
By clicking , selected Extended Property can be edited.
- (3) Delete
By clicking , selected Extended Property can be deleted.
- (4) Name
Extended Property name is indicated.
- (5) Value
Extended Property value is indicated.


The following is Extended Property procedure.

- 1 Click  and Extended Property screen is displayed.




The screenshot shows the 'User Profile' window with the 'PROPERTIES' tab selected. An 'EXTENDED PROPERTIES' dialog is open, containing a table with columns 'Name' and 'Value'. Below the table is an 'Extended Property' form with 'Name' and 'Value' input fields, each marked with a red asterisk. 'Add' and 'Cancel' buttons are at the bottom of the dialog. The 'Add' button in the dialog is highlighted with a red box.

- 2 Enter Name and Value.



The screenshot shows the 'Extended Property' dialog with the 'Name' field containing 'Age' and the 'Value' field containing '22'. Both fields are highlighted with a red box. The 'Add' button is also highlighted with a red box.


- 3 If clicking [Add], above entry is reflected into [User Profile] [PROPERTIES] screen.



The screenshot shows the 'User Profile' window with the 'PROPERTIES' tab selected. The 'EXTENDED PROPERTIES' table now contains one entry: 'Age' in the 'Name' column and '22' in the 'Value' column. This entry is highlighted with a red box. The 'Add' and 'Cancel' buttons are visible at the bottom.

The following is Extended Property edit procedure.


- 1 By clicking , selected Extended Property is displayed.



- 2 Enter Value.





- 3 By clicking [Save], above entry is reflected into [User Profile] [PROPERTIES] screen.




The following is Extended Property deletion procedure.

- 1 Select Extended Property to delete.




- 2 By clicking , selected Extended Property is deleted.



7.13 View event log (Only for administrator account)

Errors occurring in server, and so on can be referred. Click [CONFIGURATION] - [ADMINISTRATION] - [View event log]. Event log is displayed.



EVENT LEVEL	Component	Description	Date
All			
Error			
Warning			
Info			
Trace			
	DocumentMa...	Document # {0} not found.	2006/08/01 ..
	DocumentMa...	Document # {0} not found.	2006/08/01 ..
	DocumentMa...	Document # {0} not found.	2006/08/01 ..
	DocumentMa...	Document # {0} not found.	2006/08/01 ..
	ACM	null encoding specified. Assuming soap	2006/08/01 ..
	Diagnostics...	Truncated: Already Connecting for membe...	2006/08/01 ..
	ACM	null encoding specified. Assuming soap	2006/08/01 ..
	ACM	null encoding specified. Assuming soap	2006/08/01 ..

Figure7-88 View event log

Description is displayed by clicking Component.




Figure7-89 View event log

7.14 View software version

Server software version is displayed.

Click [CONFIGURATION] - [ADMINISTRATION] - [View software version].

Version is indicated as the figure below.

The figure illustrates the process of viewing the software version in the Sentinel Customer Support interface. It is divided into two parts:

Top Part: Navigation Menu

- SERVICE**
- CONFIGURATION** (highlighted)
- PERSONAL CONFIGURATION**
 - Edit your profile
- SYSTEM CONFIGURATION**
 - Manage devices
 - Manage customers
 - Manage alerts and alarms
 - Manage software packages
 - Manage connection types
 - Launch CMDB editor
- ADMINISTRATION** (highlighted)
 - View TotalAccess
 - View event log
 - View software version** (highlighted with a red box)
 - View software license

A large white arrow points from the 'View software version' option to the screenshot below.

Bottom Part: Sentinel Customer Support Screenshot

The screenshot shows the Sentinel Customer Support web application in Microsoft Internet Explorer. The user is signed in as 'hirota'. The interface includes a 'Daily Alert Summary' and 'Recent Activities' section. A 'Software Version' dialog box is open, displaying the following information:

Client Version: 5.1.106

Installed Date	Version	Description	Status
2006/07/28 12:...	5.1.102	Addresses bugs ...	Installed
2006/07/27 10:...	5.1.102	Qestra Smar...	Installed

The dialog box also shows a 'Close' button and a 'Powered by QESTRA' logo. The background interface shows a list of alerts and a sidebar with navigation options.

Figure7-90 View software version




Figure7-91 View software version

7.15 View software license (Only for administrator account)

Server software license is indicated. Click [CONFIGURATION] - [ADMINISTRATION] - [View software license]. Software license is indicated as the figure below.






Figure7-92 View software license




The image shows a 'Software License' dialog box with a table of installed licenses. The table has four columns: Component, Licensee, Serial Number, and Expire Date. A red box highlights the first eight rows of the table. A 'Close' button is located at the bottom right of the dialog box.

Component	Licensee	Serial Number	Expire Date
Control Center	Hitachi Medical ...	304	Never
Reports Server	Hitachi Medical ...	304	Never
Usage Tracking	Hitachi Medical ...	304	Never
Device Manage...	Hitachi Medical ...	304	Never
Total Access	Hitachi Medical ...	304	Never
Monitoring	Hitachi Medical ...	304	Never
Configuration M...	Hitachi Medical ...	304	Never
Software Director	Hitachi Medical ...	304	Never

Figure7-93 View software license

Chapter8 How to quit Sentinel Customer Support

- 1 Click [Sign out].



2 Click .

